



Instrumenten/ Instruments

ALIGNER PLIERS

IX887

The Torque Aligner Plier

The IX887 creates a horizontal indentation (0.75 x 1.5mm), creating a force point for individual lingual or labial root torque.



IX892

Ixion Aligner Scissors

These lightweight scissors are used for removing excess thermoformed material, and for the smooth trimming and shaping of aligners. Their streamline design, graduated curve and tapered beak makes for easy use - even in tight areas.



IX889

The Retention Aligner Plier


to create a circular indentation (1mm).



IX888

The Rotation Aligner Plier


The IX888 creates a vertical indentation / force point (0.75 x 1.5mm) on the tooth surface to optimise tooth rotation.



IX891

The Hole Punch Plier

Creates a half-circle cut-out in the aligner to give clearance for bonded buttons, tubes or brackets with hooks for elastics.



IX890

The Tear Drop Plier

Create a cut at the gingival margin, anywhere along the arch - buccally, labially or lingually - to conveniently hold elastics in place while the aligner is placed in the mouth.



LINGUAL SERIES



IX949

Slim Tweed Plier

This instrument has slimmer tips - just 0.9mm



IX950

Weingart 60° Angle

Maximum .018" (.46mm) wire.



IX952

Stepped Tweed Arch Plier

Designed for torquing and bending both square and rectangular wires.



IX953

Cinch Back

The serrated tips of this instrument hold the wire securely while bending the wire 90° distal to the buccal tube.

For use with Nickel Titanium wires
.012" - .018" (.30 x .46mm).



IX954

Cinch Back - Large Wire

Same features as ix953 but designed to work with larger wires.

For use with wires from .016" (.41mm)
up to .021" x .025" (.55 x .64mm).



IX955

Bracket Removing Plier



IX956

Mini Ligature Cutter 50° Angle Upper

Maximum Cutting Capacity, soft wire up to .012" (.30mm)



IX956s

Micro Ligature Cutter Upper 50° Angle Upper

Maximum Cutting Capacity, soft wire up to .012" (.30mm)



LINGUAL SERIES



IX911
Micro Ligature Cutter
Maximum Cutting Capacity: soft wire up to .012" (.30mm)

IX911L
Micro Long Handled Ligature Cutter
Maximum Cutting Capacity: soft wire up to .012" (.30mm)

IX958 / IX958S
Flush Cut Distal End Cutter
Standard / Small
Minimum Cutting Capacity: ix958 - .016" - ix958s .012"
Maximum Cutting Capacity: ix958 archwires up to .021" x .025"
ix958s archwires up to .016" x .025"

IX959
Long Handled Distal End Cutter
With T.C. inserts.
Minimum Cutting Capacity: .012"
Maximum Cutting Capacity: .018"x.025"

IX632
Lingual Bracket and Bondable Tube Holder

IX941
Curved Lingual Micro Mosquito

IX940
Curved Smaha
excellent for placing elastomerics to brackets


CUTTERS



IX900

Distal End Cutter

Cutting Capacity: archwires up to .021" x .025"



IX900s

Distal End Cutter

Allows you to routinely cut beta titanium, braided and co-axial wires with ease. Cutting Capacity: archwires up to .016" x .022"



IX901

Long Handled Distal End Cutter

Cutting Capacity: archwires up to .021" x .025"



IX905

Hard Wire Cutter

Minimum Cutting Capacity: .012"
Maximum Cutting Capacity: archwires up to .021" x .025"



iXion[™]
Instruments

IX902 / IX902s

Flush Cut Distal End Cutter Standard / Small

IX902


Minimum Cutting Capacity: .016"
Maximum Cutting Capacity: archwires up to .021" x .025"

IX902s

Minimum Cutting Capacity: .012"
Maximum Cutting Capacity: archwires up to .016" x .025"



Using none flush cut distal end



Using flush cut distal end



CUTTERS


IX915

Hard Wire Cutter with T.C. Inserts
Minimum Cutting Capacity: .012"
Maximum Cutting Capacity: archwires up to .021"x .025"




IX906

Hard Wire Cutter 15° Angle
Minimum Cutting Capacity: .012"
Maximum Cutting Capacity: archwires up to .021"x .025"



IX916

Hard Wire Cutter 15° Angle with T.C. Inserts
Minimum Cutting Capacity: .012"
Maximum Cutting Capacity: archwires up to .021"x .025"




IX908

Heavy Duty Hard Wire Cutter
Minimum Cutting Capacity: .014" (.36mm)
Maximum Cutting Capacity: wire up to .060" (1.5mm)




IX910

Pin Ligature Cutter 7° Angle
Maximum Cutting Capacity: soft wire up to .016"



IX911

Micro Ligature Cutter
The inserted tips ensure precision cutting of ligature wire and elastics. Maximum Cutting Capacity: soft wire up to .012"




IX911L

Micro Long Handled Ligature Cutter
Maximum Cutting Capacity: soft wire up to .012"



IX912

Standard Ligature Cutter
Maximum Cutting Capacity: soft wire up to .015"



IX913

Standard Soft Wire Ligature Cutter
Maximum Cutting Capacity: soft wire up to .012"




UTILITY PLIERS

IX800

Small Weingart


Holds round, square and rectangular wires up to .022" x .028



IX801

Small Weingart with T.C.

Inserted T.C. tips providing maximum wire control. Holds round, square and rectangular wires up to .022" x .028"



IX802

Large Weingart

Holds wire up to .030"



IX803

Large Weingart with T.C.

Holds wire up to .030"



IX804

Fine Tipped Weingart

Holds wire up to .022" x .028"



IX843

Utility Weingart with Inserted Tips

Holds wire up to .022" x .028"



IX805

How Plier (Straight)

Holds wire up to .022" x .028"



IX806

How Plier (Angled)

Holds wire up to .022" x .028"




IX820

Crimpable Archwire Hook Plier



IX834

Crimping Plier



IX825

Posterior Band Remover (Short)

Replacement Tips 3/16" (5mm) - IX895S

Replacement Tips 1/8" (4mm) - IX896S



UTILITY PLIERS

IX826

Posterior Band Remover (Long)
Replacement Tips 3/16" (5mm) - IX895S
Replacement Tips 1/8" (4mm) - IX896S



IX838

Posterior Band Remover Plier
Replacement Tips. Pack of 10.



IX827

Direct Bond Remover (Straight)



IX828

Direct Bond Remover (Angled 60°)




IX829

Direct Bond Remover



IX840

Small Adams Plier
Effective for wires up to .030"



IX841

Force Module Separator



IX895S

Replacement Tips 3/16" (5mm)
for IX825/826, Pk 3



IX896S

Replacement Tips 1/8" (4mm)
for IX825/826, Pk 3



IX838S

Posterior band remover tips
for IX838 Posterior Band Remover




WIRE BENDING PLIERS

IX807

Light Wire Plier

Tips bend archwires from .012" to .020".



IX809

Light Wire Plier with Cutter

Tips bend archwires from .012" to .020".
Cutter at the base of beaks cuts archwires from .014" to .020".



IX810


Bird Beak Plier



IX812

Tweed Loop Forming Plier with Cutter

The tip has three sections of 1.2mm, 1.8mm and 2.0mm.
For use with round and square wire up to .021 x .025".
The built in cutter allows wire to be cut from .014" to .020".



IX813

Tweed Arch Forming Plier

For use with wires up to .021" x .025".



IX821

Long Beak Light Wire Plier


For use with round and square wire up to .020" (.51mm).



WIRE BENDING PLIERS

IX872

Bird Beak Plier with T.C.



IX873

Torquing Plier Set



IX819

Small 3-Jaw Plier
For working with wire up to .020".



IX814

O'Brien Loop Forming Plier
The perfect plier for making Accurate, Omega and Closing loops. Can be used with round & square wire up to .021" x .025" (.55 x .64mm).



IX815

Nance Plier
Precision loop forming with four step tips 3mm, 4mm, 5mm and 6mm. Excellent for round and square wires up to .022".




IX816

3-Jaw Plier
For use with wire up to .030".



IX817

De La Rosa
For use with wires up to .021" x .025".



WIRE BENDING PLIERS



IX830

Hollow Chop



IX975

Intra Oral Detailing Plier 0.75mm Step

For use up to max .018" x.018" wire.



IX823

Fine 3-Jaw Plier

Can be used on Nickel Titanium and
Stainless Steel wire up to 0.20".



IX976

Intra Oral Detailing Plier 0.5mm Step

For use up to max
.018" x.018" wire.



HAND INSTRUMENTS & HEMOSTATS



IX650
Probe Single Ended



IX651
Probe Double Ended



IX652
Ligature Tucker



IX653
Mitchell Trimmer



IX656
Cement Spatula



IX657
Flat Plastic



IX654
Mirror Handle



IX655
Mirror Heads Pk 12

HAND INSTRUMENTS & HEMOSTATS

IX850

Direct Bond Bracket Holder



IX630

Direct Bond Holder with Positioning Blade



IX631

Deluxe Direct Bond Holder



IX632

Lingual Bracket and Bondable Tube Placer



IX633

Curved Bracket Holder



IX845

Mershon Band Pusher



IX861

Mathieu Needle Holder with T.C.
Designed for twisting
ligature wire.



IX862

Mathieu Needle Holder
Designed for twisting
ligature wire.



IX863


Smaha (Straight)
Fine tips excellent for placing
elastomeric ligatures.



HAND INSTRUMENTS & HEMOSTATS

IX855

Mosquito Forceps



IX860

Fine Tip Mathieu Needle Holder



IX940

Curved Smaha


Excellent for placing elastomerics to lingual brackets.



IX941

Curved Lingual Micro Mosquito

Used to place elastic ligatures.



IX866

Bracket Height Gauge

3.5mm-5mm .018

For precise bracket positioning of .018 brackets.



IX867

Bracket Height Gauge

3.5mm-5mm .022

For precise bracket positioning of .022 brackets.



HAND INSTRUMENTS & HEMOSTATS

IX892

Ixion Aligner Scissors



IX101

Small Rack
Autoclavable
Stainless Steel Rack
size L 110mm x H 120mm



IX103

Large Rack Autoclavable Stainless Steel Rack
size L 310mm x H 125mm



IX102

Medium Rack
Autoclavable
Stainless Steel Rack
size L 155mm x H 120mm



DB04-000/1

Clear Upright Instrument Rack
Keep your instruments organised
will hold up to 15 instruments



LABORATORY PLIERS

IX164

Adams Plier

Bends hard wire up to .039" (1.00mm) max.



IX165

Spring Former Plier

With tungsten Carbide inserted tips.
Bends hard wire up to 0.28" (0.7mm) max.



IX700

Universal General Use Plier

Both Tips are serrated. Bends hardwire up to .039" (1.0mm) max. Cuts hardwire up to .028" (0.7mm) max.



IX702

Three Jaw Plier

For hard wire up to .024" (1mm) max.



IX703

Young Plier


For labial arches and canine loops.
For hard wire up to .028" (0.7mm) max.



IX705

Lingual Arch Bending Plier

Grooved to hold all larger wire sizes, jaw width is .120" (3.05mm). Also ideal for double and triple back bends. Bends round, square and rectangular wires up to .036" (0.91mm) max.



IX708

Optical Plier

For hardwire up to .028" (0.7mm) max.



IX711

Heavy Duty Bird Beak

Bending of wire up to .040" 1mm.



INSTRUMENTS

BRACKET DELUXE HOLDER

DB05-0421



STANDARD BRACKET HOLDER

DB05-0420



BONDABLE TUBE HOLDER

DB05-0423



CEMENT SPATULA

Multifunctional orthodontic hand instrument with precision manufactured tip.

DB05-0325



COLLEGE TWEEZERS

DB05-0319



PEDO CURVED SCISSORS

DB05-0331



OVERBITE RULER

Aluminium ruler for easy measuring of overbite.


DB05-0001



ARCH FORMER - Tweed Style

Will bend wires up to .022" x .028".

DB05-0712



DOUBLE ENDED PROBE

DB05-0312



INSTRUMENTS


SPECTRUM STEP PLIER

800-276 (1mm)
800-277 (3/4mm)
800-278 (1/2mm)




SWIVEL HEAD BRACKET Height Gauge

800-068



SPECTRUM TORQUING PLIER SET w/Wide Key

800-130
800-131
replacement key



LARGE MOUTH MIRROR

800-050 38,10mm



BOONE BRACKET Gauge

800-034



DOUBLE ENDED ELASTIC REMOVER

800-069



BAND SEATER/BITE STICK


autoclavable and cold sterilization

800-0182



REVERSE ACTION KEAT MOLAR TUBE PLACER

800-039L



CINCH BACK / DISTAL BENDER

800-075



BRACKET POSITIONER

800-046



5-IN-1 BRACKET PLACER & ALIGNER

800-041



INSTRUMENTS

SHORTY TWIST HOLDER

800-133



LITE DOUBLE ENDED LIGATURE DIRECTOR

800-062 (1/4")



LITE BAND PUSHER / SCALER

(3/8") 800-064W



LITE LIGATURE REMOVER / DIRECTOR

(1/4") 800-065 / (3/8") 800-065W



MATHIEU NEEDLE HOLDER w/Safety Latch-Regular Tip

A wider serrated tip to hold ligature wire and Kobyashi Hooks firmly during placement.

800-122



MATHIEU NEEDLE HOLDER w/Safety Latch-Fine Tip

A finer serrated tip to securely place elastic ligature and wires.

800-120



MATHIEU NEEDLE HOLDER w/Safety Latch-Ultra Fine Tip

This ultra-fine tip is designed specifically to place elastic ligatures securely and more easily by taking up less of the donut.

800-121



INSTRUMENTS

BRACKETIZER Blue Handle

800-076



50-2030-0000

For Experience metal/rhodium



50-0032-0000

For Experience Mini metal/rhodium



DB05-0311

Double ended tucker



LM411-415ES

LM Ligature Director



50-2130-H000

The universal Experience instrument's (handle only)

50-2130-B000

Additional blades
(Verpakt per 2 / Emballé par 2)



50-0030-0000

For Experience metal/
ceramic and rhodium

HORIZONTAL PLIER RACK w/ Cover (12 Pliers)

700-074



UPRIGHT PLIER RACK w/Cover (12 pliers)

700-071



OIL PEN van Carl Martin

990MF20



PLIER CAROUSEL w/ Center Tube

700-075



SILICONE COLOR CODED INSTRUMENT RINGS

(80)

700-050



