



## Elastieken/ Élastiques



## ORCHESTRA INTRA & EXTRAORAL ELASTICS

Force/Size	EXTRA-ORAL			
	2.5 oz (70,88 g)	4.5 oz (113,40 g)	6.5 oz (170,10 g)	8,5 oz (240,98g)
3 mm (1/8")	85-0321-2000 <b>Harp</b>	85-0331-2000 <b>Xylophone</b>	85-0341-2000 <b>Tuba</b>	-
4 mm (3/16")	85-0421-2000 <b>Tambourine</b>	85-0431-2000 <b>Piano</b>	85-0441-2000 <b>Trumpet</b>	85-0461-2000 <b>Trombone</b>
4 mm (3/16") non latex	-	-	85-0440-0000 <b>Rock-n-Roll</b>	-
6 mm (1/4")	85-0621-2000 <b>Oboe</b>	85-0631-2000 <b>Flute</b>	85-0641-2000 <b>Baritone</b>	85-0661-2000 <b>French Horn</b>
6 mm (1/4") non latex	-	85-0630-0000 <b>Hip-Hop</b>	85-0640-0000 <b>Classical</b>	-
8 mm (5/16")	85-0821-2000 <b>Cello</b>	85-0831-2000 <b>Clarinet</b>	85-0841-2000 <b>Saxophone</b>	85-0861-2000 <b>Guitar</b>
10 mm (3/8")	85-1021-2000 <b>Cymbals</b>	85-1031-2000 <b>Tympani</b>	-	85-1061-2000 <b>Violin</b>

## EXTREME ELASTICS



Force/Size	Light 3.5 oz	Medium 4.5 oz	Heavy 6.0 oz
3,2 mm (1/8")	M931-01 <b>Kayaking</b>	M931-06 <b>Surfing</b>	M931-11 <b>Jet Skiing</b>
4,8mm (3/16")	M931-02 <b>Motorcross</b>	M931-07 <b>Skiing</b>	M931-12 <b>Bungeejumping</b>
6,4mm (1/4")	M931-03 <b>Rock Climbing</b>	M931-08 <b>Skydiving</b>	M931-13 <b>Snowmobiling</b>
7,9mm (5/16")	M931-04 <b>Skateboarding</b>	M931-09 <b>BMX Freestyle</b>	M931-14 <b>Snowboarding</b>
9,5mm (3/8")	M931-05 <b>Water Skiing</b>	M931-10 <b>Windsurfing</b>	M931-15 <b>Wakeboarding</b>

### ELASTIC ATTACHER



100/pk  
Assorted colors **800-026**

### ELASTIC STORAGE KEY RINGS



50/pkg.  
Neon assorted **ESC050**  
Basic assorted **ESC150**

# COLOR CHART for elastomeric Ties & Chain

CLEAR	GLOW	PEARL	METALLIC RED	RUBY RED	ORANGE
WHITE	GLOW BLUE	PEARL GREEN	METALLIC GOLD	MAROON	FLUORESCENT ORANGE
BLACK	TOOTH	PEARL BLUE	METALLIC PURPLE	LAVENDER	FUCHSIA
GRAY	NEON GREEN	GREEN	TRANSPARENT GREEN	TURQUOISE	PINK
SMOKE	DEEP SKY BLUE	NAVY BLUE	METALLIC BLUE	ROYAL BLUE	FLUORESCENT PINK
METALLIC SILVER	LIGHT SKY BLUE	YELLOW			




## LIGITIEBOOK


- 20 ligaturen per stick
- 50 sticks verpakking
- 1000 ligaturen per verpakking

- 20 ligatures par bâtonnet
- paquet de 50 bâtonnets
- 1000 ligatures par paquet

## DYNA-LINK™ ELASTOMERIC CHAIN Latex-Free



- 1 rol chain 4,5m, .0255" dikte
- 1 rouleau de chaînette 4,5m, .0255" d'épaisseur


		Select-Tie	Dyna-Link™ Chain		
					
			Continuous	Short	Long
	Pearl	OCSLPL	CDLPL	SDLPL	LDLPL
	Glow	OCSLGL	CDLGL	SDLGL	LDLGL
	Neon Green	OCSLNG	CDLNG	SDLNG	LDLNG
	Pearl Green	OCSLPG	CDLPG	SDLPG	LDLPG
	Glow Blue	OCSLGB	CDLGB	SDLGB	LDLGB
	Pearl Blue	OCSLPB	CDLPB	SDLPB	LDLPB
	Lavender	OCSLLV	CDLLV	SDLLV	LDLLV
	Maroon	OCSLMA	CDLMA	SDLMA	LDLMA
	Clear	OCSLCL	CDLCL1	SDLCL1	LDLCL1
	Gray	OCSLGY	CDLGY	SDLGY	LDLGY
	Metallic Silver	OCSLMS	CDLSV	SDLSV	LDLSV
	Metallic Gold	OCSLGD	CDLGD	SDLGD	LDLGD
	Metallic Purple	OCSLPP	CDLPP	SDLPP	LDLPP
	Metallic Red	OCSLMR	CDLMR	SDLMR	LDLMR
	Metallic Blue	OCSLMB	CDLMB	SDLMB	LDLMB
	Royal Blue	OCSLRB	CDLRB	SDLRB	LDLRB
	Lt. Sky Blue	OCSLLS	CDLLS	SDLLS	LDLLS
	Deep Sky Blue	OCSLDS	CDLDS	SDLDS	LDLDS
	Navy Blue	OCSLNB	CDLNB	SDLNB	LDLNB
	Green	OCSLGR	CDLGR	SDLGR	LDLGR
	Transp. Green	OCSLTG	CDLTG	SDLTG	LDLTG
	Ruby Red	OCSLRR	CDLRR	SDLRR	LDLRR
	Turquoise	OCSLTQ	CDLTQ	SDLTQ	LDLTQ
	White	OCSLWT	CDLWT	SDLWT	LDLWT
	Fuchsia	OCSLFU	CDLFU	SDLFU	LDLFU
	Black	OCSLBK	CDLBK	SDLBK	LDLBK
	Orange	OCSLOR	CDLOR	SDLOR	LDLOR
	Fluor. Orange	OCSLFO	CDLFO	SDLFO	LDLFO
	Yellow	OCSLYW	CDLYW	SDLYW	LDLYW
	Pink	OCSLPK	CDLPK	SDLPK	LDLPK
	Fluor. Pink	OCSLFP	CDLFP	SDLFP	LDLFP
	Tooth	OCSLTC	CDLTC	SDLTC	LDLTC
	Smoke	OCSLSM	CDLSM	SDLSM	LDLSM

# ELASTICS LIGATURES



colour  
boutique

	Elastic Ligatures	Closed Chain	Open Chain	Open Long Chain
● Red	DB03-3000	DB03-3100	DB03-3200	DB03-3300
● Orange	DB03-3001	DB03-3101	DB03-3201	DB03-3301
● Silver	DB03-3002	DB03-3102	DB03-3202	DB03-3302
● Neon Pink	DB03-3003	DB03-3103	DB03-3203	DB03-3303
● Pink	DB03-3004	-	-	-
● Baby Pink	DB03-3005	DB03-3105	DB03-3205	DB03-3305
● Deep Pink	DB03-3006	DB03-3106	DB03-3206	DB03-3306
● Lavender	DB03-3007	DB03-3107	DB03-3207	DB03-3307
● Neon Green	DB03-3008	DB03-3108	DB03-3208	DB03-3308
● Purple	DB03-3009	DB03-3109	DB03-3209	DB03-3309
● Baby Blue	DB03-3010	-	-	-
● Black	DB03-3011	DB03-3111	DB03-3211	DB03-3311
● Metallic Blue	DB03-3012	DB03-3112	DB03-3212	DB03-3312
● Navy Blue	DB03-3013	DB03-3113	DB03-3213	DB03-3313
● Blue	DB03-3014	DB03-3114	DB03-3214	DB03-3314
● Metallic Purple	DB03-3015	-	-	-
● White	DB03-3016	DB03-3116	DB03-3216	DB03-3316
● Clear	DB03-3017	DB03-3117	DB03-3217	DB03-3317
● Cyan	DB03-3018	DB03-3118	DB03-3218	DB03-3318
● Turquoise	DB03-3019	DB03-3119	DB03-3219	DB03-3319
● Lilac	DB03-3020	DB03-3120	DB03-3220	DB03-3320
● Green	DB03-3021	-	-	-
● Racing Green	DB03-3022	-	-	-
● Aubergine	DB03-3023	-	-	-
● Sky Blue	DB03-3024	-	-	-
● Ice White	DB03-3025	-	-	-
● Mixed pack	DB03-3026	-	-	-
● Yellow	DB03-3027	DB03-3121	DB03-3221	DB03-3321
● Gold	DB03-3028	DB03-3122	DB03-3222	DB03-3322



## Elastic ligatures

- 26 ligatures per stick
- 40 sticks verpakking
- 1040 ligaturen per verpakking

- 26 ligatures par bâtonnet
- paquet de 40 bâtonnets
- 1040 ligatures par paquet

## Chain

- 1 rol chain 4,5m

- 1 rouleau de chaînette 4,5m



## GC ORTHOCHAIN

	Short - 2m	Medium - 2m	Long - 2m
○ Clear	86-1010-000	86-1020-000	86-1030-000
● Grey	86-2010-000	86-2020-000	86-2030-000

# ELASTICS


## ELASTOMERIC LIGATURE KIT

- Lege kit / Boîte vide: **DB03-L/box**



## CHAIN DISPENSER

- 18 spools **700-065**



## TIES/CHAIN COMBO DISPENSER **700-062**



## ELASTOMERIC SEPARATORS

**325-302**



## ELASTOMERIC THREAD & TUBING

- Tube Gray  
25 ft. Spool - .025 - **340-210**
- Tube Clear  
25 ft. Spool - .025 - **340-205**
- Thread Gray  
25 ft. Spool - .025 - **340-110**
- Thread Clear  
25 ft. Spool - .025 - **340-105**
- Slip-not Solid - .025  
**SNTGY025** grey / **SNTCL025** clear



## ROTATION WEDGES


100/pkg

- Clear **270-021**
- Grey **270-022**



## DENTALASTICS® SEPARATORS

- 1000 pcs - BLUE - **774-200-01**




## ARCH SLEEVE PROTECTIVE TUBING

**Latex-Free** - 10 ft. Per spool - .027" - .027"x.027"

- Grey **340-310**
- Clear **340-305**



## SLINGTIES™




Verpakking / Emballage  
- 104 x 6mm Slingties  
- 104 x 9mm Slingties™  
**300-000**



De Slingties sluiten de diastemen en verankeren 2 tanden.  
Les Slingties ferment le diastème et forme un ancrage avec 2 dents.



De Slingties maken ruimte en brengt de tand naar de draad.  
Les Slingties créent de l'espace et amènent la dent vers le fil.



De Slingtie trekt de hoektaand naar de draad.  
Le Slingtie tire la canine vers le fil.

