



**Kleef- & Weldable tubes/  
Tubes à coller & à souder**  
**TUBES**



# SELF-LIGATING Tubes



## Experience Self-ligating tubes

### MAXILLARY

Tooth 6 & 7	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
Experience SL tube, Roth	-10°	0°	14°D	20-1160-0010	20-1260-0010	20-2160-0010	20-2260-0010
Experience SL tube, Small Base Roth	-10°	0°	14°D	-	-	20-2160-S010	20-2260-S010
Experience SL tube, MBT	-14°	0°	10°D	20-1160-0020	20-1260-0020	20-2160-0020	20-2260-0020
Experience SL tube, Small Base MBT	-14°	0°	10°D	20-1160-S020	20-1260-S020	20-2160-S020	20-2260-S020
Experience double SL tube, Extra Torque	-20°	0°	20°D	20-1160-0030	20-1260-0030	20-2160-0030	20-2260-0030

### MANDIBULAR

Tooth 6 & 7	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
Experience SL tube, Roth	-25°	0°	4°D	20-1460-0010	20-1360-0010	20-2460-0010	20-2360-0010
Experience SL tube, Small Base Roth	-25°	0°	4°D	-	-	20-2460-S010	20-2360-S010
Experience SL tube, MBT	-20°	0°	0°D	20-1460-0020	20-1360-0020	20-2460-0020	20-2360-0020
Experience SL tube, Small Base MBT	-20°	0°	0°D	20-1460-S020	20-1360-S020	20-2460-S020	20-2360-S020
Experience SL tube, Extra Torque	-14°	0°	4°D	-	-	20-2460-0032	20-2360-0032

### MANDIBULAR

Tooth 6 & 7	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
Experience DOUBLE SL tube, Extra Torque	-10°	0°	4°D	20-1460-0030	20-1360-0030	20-2460-0030	20-2360-0030

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# LEGEND LP Tubes

## Conventional tubes

Tooth	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
U6 , Roth	-14°	0°	14°	40-1160-0114	40-1260-0114	40-2160-0114	40-2260-0114
U7, Roth	-14°	0°	14°	40-1170-0114	40-1270-0114	40-2170-0114	40-2270-0114
L6, Roth	-25°	0°	4°	40-1460-0114	40-1360-0114	40-2460-0114	40-2360-0114
L7, Roth	-25°	0°	6°	40-1470-0114	40-1370-0114	40-2470-0114	40-2370-0114
U6 , MBT	-14°	0°	10°	40-1160-0106	40-1260-0106	40-2160-0106	40-2260-0106
U6 , MBT	-14°	0°	0°	40-1160-0107	40-1260-0107	40-2160-0107	40-2260-0107
U7, MBT	-14°	0°	10°	40-1170-0106	40-1270-0106	40-2170-0106	40-2270-0106
U7, MBT	-14°	0°	0°	40-1170-0107	40-1270-0107	40-2170-0107	40-2270-0107
L6, MBT	-20°	0°	0°	40-1460-0106	40-1360-0106	40-2460-0106	40-2360-0106
L7, MBT	-10°	0°	0°	40-1470-0106	40-1370-0106	40-2470-0106	40-2370-0106



# LEGEND LE Tubes



## LEGEND LE Tubes (Large Entrance)

### Second Molar Bondable

Non-convertible single tubes for second molars, the funnel entry of these tubes facilitates access to second molars and its anatomical base fits the teeth perfectly making positioning easy.

### MAXILLARY

#### Type Roth

Dimensions	Torque	Rotation	P/N	
			R	L
.018"x.025"	-14°	14°D	40-1170-0099	40-1270-0099
.022"x.028"	-14°	14°D	40-2170-0099	40-2270-0099

#### Type Andrews

Dimensions	Torque	Rotation	P/N	
			R	L
.018"x.025"	-9°	7°D	40-1170-0100	40-1270-0100
.022"x.028"	-9°	7°D	40-2170-0100	40-2270-0100

### MANDIBULAR

#### Type Roth

Dimensions	Torque	Rotation	P/N	
			R	L
.018"x.025"	-30°	6°D	40-1470-0099	40-1370-0099
.022"x.028"	-30°	6°D	40-2470-0099	40-2370-0099

#### Type Andrews

Dimensions	Torque	Rotation	P/N	
			R	L
.018"x.025"	-25°	0°	40-1470-0100	40-1370-0100
.022"x.028"	-25°	0°	40-2470-0100	40-2370-0100

#### Type MBT\*

Dimensions	Torque	Rotation	P/N	
			R	L
.018"x.025"	-10°	0°	40-1470-0101	40-1370-0101
.022"x.028"	-10°	0°	40-2470-0101	40-2370-0101

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# MINI Tubes




## MINI ERUPTION TUBE

for 2<sup>nd</sup> Molar, low profile bondable

Dimensions	Torque	Rotation	Length mm	P/N	
				R	L
.018"x.025"	-10°	0°D	2.0	40-1170-0122	40-1270-0122
.022"x.028"	-10°	0°D	2.0	40-2170-0122	40-2270-0122
.018"x.025"	-25°	0°D	2.0	40-1470-0122	40-1370-0122
.022"x.028"	-25°	0°D	2.0	40-2470-0122	40-2370-0122

# AXCESS Tubes



## AXCESS Tubes

These entry point buccal tubes are manufactured in one-piece using MIM (Metal Injection Molded) technology. AXCESS tubes are available in bondable and weldable versions.

### MAXILLARY

Double Non-Convertible HG .045" 1<sup>st</sup> Molar

Dimensions	Torque	Rotation	R/L	P/N	
				Weldable	Bondable
.018"x.025"	-10°	6°D	R	41-1160-0023	—
			L	41-1260-0023	—
.022"x.028"	-10°	6°D	R	41-2160-0023	—
			L	41-2260-0023	—
.018"x.025"	0°	6°D	R	41-1160-0024	—
			L	41-1260-0024	—
.022"x.028"	0°	6°D	R	41-2160-0024	—
			L	41-2260-0024	—

Single Non-Convertible 1<sup>st</sup> Molar

Dimensions	Torque	Rotation	R/L	P/N	
				Weldable	Bondable
.022"x.028"	-10°	6°D	R	41-2160-0025	40-2160-0025
			L	41-2260-0025	40-2260-0025
.018"x.025"	0°	6°D	R	41-1160-0026	40-1160-0026
			L	41-1260-0026	40-1260-0026
.022"x.028"	0°	6°D	R	—	40-2160-0026
			L	—	40-2260-0026

Single Non-Convertible 2<sup>nd</sup> Molar

Dimensions	Torque	Rotation	R/L	P/N	
				Weldable	Bondable
.018"x.025"	-10°	6°D	R	—	40-1170-0088
			L	—	40-1270-0088
.018"x.025"	0°	0°	R	—	40-1170-0090
			L	—	40-1270-0090
.018"x.025"	-10°	0°	R	—	40-1170-0091
			L	—	40-1270-0091

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# AXCESS Tubes



## AXCESS Tubes

### Single Non-Convertible 2<sup>nd</sup> Molar

Dimensions	Torque	Rotation	R/L	P/N	
				Weldable	Bondable
.022"x.028"	-10°	6°D	R	41-2170-0025	40-2170-0025
			L	41-2270-0025	40-2270-0025
.022"x.028"	-10°	6°D	R	—	40-2170-0088
			L	—	40-2270-0088
.022"x.028"	0°	0°	R	—	40-2170-0090
			L	—	40-2270-0090
.022"x.028"	-10°	0°	R	—	40-2170-0091
			L	—	40-2270-0091

## MANDIBULAR

### Single Non-Convertible 1<sup>st</sup> Molar

Dimensions	Torque	Rotation	R/L	P/N	
				Weldable	Bondable
.022"x.028"	-25°	6°D	R	41-2460-0025	40-2460-0025
			L	41-2360-0025	40-2360-0025
.018"x.025"	-25°	6°D	R	41-1460-0025	-
			L	41-1360-0025	-
.022"x.028"	0°	6°D	R	41-2460-0026	40-2460-0026
			L	41-2360-0026	40-2360-0026
.018"x.025"	0°	6°D	R	41-1460-0026	40-1460-0026
			L	41-1360-0026	40-1360-0026

### Single Non-Convertible 2<sup>nd</sup> Molar

Dimensions	Torque	Rotation	R/L	P/N	
				Weldable	Bondable
.018"x.025"	0°	0°	R	—	40-1470-0090
			L	—	40-1370-0090
.018"x.025"	-25°	0°	R	—	40-1470-0091
			L	—	40-1370-0091
.022"x.028"	-25°	6°D	R	41-2470-0025	40-2470-0025
			L	41-2370-0025	40-2370-0025
.022"x.028"	-25°	6°D	R	—	40-2470-0089
			L	—	40-2370-0089
.022"x.028"	0°	0°	R	—	40-2470-0090
			L	—	40-2370-0090
.022"x.028"	-25°	6°D	R	—	40-2470-0091
			L	—	40-2370-0091

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Maxillary Single Tubes Convertible



Type Roth, 1<sup>st</sup> Molar

Dimensions	Torque	Rotation	Length mm	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-14°	14°D	4.0	R	41-1160-0061	40-1160-0061
				L	41-1260-0061	40-1260-0061
.022"x.028"	-14°	14°D	4.0	R	41-2160-0061	40-2160-0061
				L	41-2260-0061	40-2260-0061



Type Roth, 1<sup>st</sup> Molar

Dimensions	Torque	Rotation	Length mm	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-10°	14°D	4.0	R	41-1160-0012	40-1160-0012
				L	41-1260-0012	40-1260-0012
.022"x.028"	-10°	14°D	4.0	R	41-2160-0012	40-2160-0012
				L	41-2260-0012	40-2260-0012



Type MBT\*, 1<sup>st</sup> Molar

Dimensions	Torque	Rotation	Length mm	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-14°	10°D	3.8	R	41-1160-0022	40-1160-0022
				L	41-1260-0022	40-1260-0022
.022"x.028"	-14°	10°D	3.8	R	41-2160-0022	40-2160-0022
				L	41-2260-0022	40-2260-0022



Type MBT\*, 1<sup>st</sup> Molar

Dimensions	Torque	Rotation	Length mm	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-14°	10°D	4.2	R	41-1160-0110	40-1160-0110
				L	41-1260-0110	40-1260-0110
.022"x.028"	-14°	10°D	4.2	R	41-2160-0110	40-2160-0110
				L	41-2260-0110	40-2260-0110



Type Andrews, 1<sup>st</sup> Molar

Dimensions	Torque	Rotation	Length mm	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-10°	7°D	3.6	R	41-1160-0034	40-1160-0034
				L	41-1260-0034	40-1260-0034
.022"x.028"	-10°	7°D	3.6	R	41-2160-0034	40-2160-0034
				L	41-2260-0034	40-2260-0034

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Maxillary Single Tubes Non-Convertible



GET 1<sup>st</sup> Molar



Dimensions	Torque	Rotation	Length mm	R/L	P/N	
					Weldable	Bondable
.022"x.028"	-20°	10°D	4.0	R	—	40-2160-0510
				L	—	40-2260-0510
.022"x.028"	-20°	0°	4.0	R	—	40-2160-0511
				L	—	40-2260-0511



GET 2<sup>nd</sup> Molar



Dimensions	Torque	Rotation	Length mm	R/L	P/N	
					Weldable	Bondable
.022"x.028"	-20°	0°	4.0	R	—	40-2170-0120
				L	—	40-2270-0120



Type Roth Distal Extension, 2<sup>nd</sup> Molar

Dimensions	Torque	Rotation	Length mm	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-14°	14°D	4.0	R	—	40-1170-0013
				L	—	40-1270-0013
.022"x.028"	-14°	14°D	4.0	R	—	40-2170-0013
				L	—	40-2270-0013



Type Roth Distal Extension, 2<sup>nd</sup> Molar

Dimensions	Torque	Rotation	Length mm	R/L	P/N	
					Weldable	Bondable
.022"x.028"	-10°	14°D	3.0	R	41-2170-0039	40-2170-0039
				L	41-2270-0039	40-2270-0039



Type Roth, 2<sup>nd</sup> Molar

Dimensions	Torque	Rotation	Length mm	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-14°	14°D	4.2	R	—	40-1170-0014
				L	—	40-1270-0014
.022"x.028"	-14°	14°D	4.2	R	41-2170-0014	40-2170-0014
				L	41-2270-0014	40-2270-0014

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs




# Maxillary Single Tubes Non-Convertible



Type Roth, 2<sup>nd</sup> Molar

						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-14°	0°	4.0	R	—	40-1170-0063	
				L	—	40-1270-0063	
.022"x.028"	-14°	0°	4.0	R	41-2170-0063	40-2170-0063	
				L	41-2270-0063	40-2270-0063	



Type Roth, 2<sup>nd</sup> Molar

						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-14°	14°D	4.0	R	41-1170-0012	40-1170-0012	
				L	41-1270-0012	40-1270-0012	
.022"x.028"	-14°	14°D	4.0	R	41-2170-0012	-	
				L	41-2270-0012	-	



Type MBT\*, 2<sup>nd</sup> Molar

						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-14°	10°D	4.2	R	41-1170-0022	40-1170-0022	
				L	41-1270-0022	40-1270-0022	
.022"x.028"	-14°	10°D	4.2	R	41-2170-0022	40-2170-0022	
				L	41-2270-0022	40-2270-0022	



Type Alexander Distal Extension, 2<sup>nd</sup> Molar

						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Weldable	Bondable	
.022"x.028"	-10°	7°D	4.0	R	41-2170-0075	—	
				L	41-2270-0075	—	



Type Andrews, 2<sup>nd</sup> Molar


						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-10°	7°D	4.0	R	41-1170-0034	40-1170-0034	
				L	41-1270-0034	40-1270-0034	
.022"x.028"	-10°	7°D	4.0	R	41-2170-0034	40-2170-0034	
				L	41-2270-0034	40-2270-0034	

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



### Type Andrews, 2<sup>nd</sup> Molar

						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-10°	7°D	3.0	R	41-1170-0048	40-1170-0048	
				L	41-1270-0048	40-1270-0048	
.022"x.028"	-10°	7°D	3.0	R	41-2170-0048	40-2170-0048	
				L	41-2270-0048	40-2270-0048	



### Type Andrews Wide Base, 2<sup>nd</sup> Molar

						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-10°	7°D	3.0	R	—	40-1170-0049	
				L	—	40-1270-0049	
.022"x.028"	-10°	7°D	3.0	R	—	40-2170-0049	
				L	—	40-2270-0049	




### Type Andrews Wide Entry, 2<sup>nd</sup> Molar

						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-10°	7°D	3.5	R	41-1170-0018	40-1170-0018	
				L	41-1270-0018	40-1270-0018	
.022"x.028"	-10°	7°D	3.5	R	41-2170-0018	40-2170-0018	
				L	41-2270-0018	40-2270-0018	



### Type Planché, 2<sup>nd</sup> Molar

						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Weldable	Bondable	
.022"x.028" -	20°	14°D	4.2	R	41-2170-0071	—	
				L	41-2270-0071	—	



### Type Polar, 2<sup>nd</sup> Molar

						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-25°	6°D	4.0	R	—	40-1170-0029	
				L	—	40-1270-0029	



### Mini Eruption Tube, 2<sup>nd</sup> Molar

						P/N	
Dimensions	Torque	Rotation	Length mm	R/L	Without hook	With hook	
.018"x.025"	-10°	0°D	2.0	R	40-1170-0112	40-1170-0019	
				L	40-1270-0112	40-1270-0019	
.022"x.028"	-10°	0°D	2.0	R	40-2170-0112	40-2170-0019	
				L	40-2270-0112	40-2270-0019	

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Maxillary Double Tubes 1<sup>st</sup> Molar Convertible



Type Roth Convertible No HG

Dimensions	Torque	Rotation Slot & Aux.	R/L	P/N	
				Weldable	Bondable
.018"x.025" / .018"x.025"	-14°	14°D	R	41-1160-0011	40-1160-0011
			L	41-1260-0011	40-1260-0011
.022"x.028" / .018"x.025"	-14°	14°D	R	41-2160-0011	40-2160-0011
			L	41-2260-0011	40-2260-0011




Type Roth Convertible No HG

Dimensions	Torque	Rotation Slot	R/L	P/N	
				Weldable	Bondable
.018"x.025" / .018"x.025"	-14°	14°D	R	41-1160-0055	40-1160-0055
			L	41-1260-0055	40-1260-0055



Type MBT\* Convertible No HG

Dimensions	Torque	Rotation Slot	R/L	P/N	
				Weldable	Bondable
.018"x.025"	-14°	10°D	R	41-1160-0104	40-1160-0104
			L	41-1260-0104	40-1260-0104
.022"x.028"	-14°	10°D	R	41-2160-0104	40-2160-0104
			L	41-2260-0104	40-2260-0104



Type MBT\* Convertible No HG

Dimensions	Torque	Rotation Slot	R/L	P/N	
				Weldable	Bondable
.018"x.025"	-14°	0°D	R	41-1160-0105	40-1160-0105
			L	41-1260-0105	40-1260-0105
.022"x.028"	-14°	0°D	R	41-2160-0105	40-2160-0105
			L	41-2260-0105	40-2260-0105

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



#### Extra Torque Convertible No HG

Dimensions	Torque	Rotation Slot	R/L	P/N	
				Weldable	Bondable
.022"x.028"	-20°	20°D	R	41-2160-0030	-
			L	41-2260-0030	-



#### Type Roth Convertible - HG .045"/Gingival with HG

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-14°	14°D	14°D	R	41-1160-0054	—
				L	41-1260-0054	—
.022"x.028"	-14°	14°D	14°D	R	41-2160-0054	—
				L	41-2260-0054	—



#### Type MBT\* Convertible - HG .045"/Occlusal With HG

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-14°	10°D	0°	R	41-1160-0021	—
				L	41-1260-0021	—
.022"x.028"	-14°	10°D	0°	R	41-2160-0021	—
				L	41-2260-0021	—



#### Type MBT\* Convertible - HG .045"/Gingival With HG

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.022"x.028"	-14°	10°D	10°D	R	41-2160-0020	-
				L	41-2260-0020	-

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Maxillary Double Tubes 1<sup>st</sup> Molar



## Type BiDimensional HG Convertible .051"/Occlusal

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.022"x.028"	0°	6°D	0°	R	41-2160-0080	—
				L	41-2260-0080	—



## Type BiDimensional HG Convertible .051"/Occlusal

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.022"x.028"	0°	14°D	7°D	R	41-2160-0083	—
				L	41-2260-0083	—



## Type Alexander Convertible - HG .045"/Occlusal

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-10°	14°D	7°D	R	41-1160-0073	—
				L	41-1260-0073	—
.022"x.028"	-10°	14°D	7°D	R	41-2160-0073	—
				L	41-2260-0073	—



## Type Andrews Convertible - HG .045"/Occlusal

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-10°	7°	7°D	R	41-1160-0079	40-1160-0079
				L	41-1260-0079	40-1260-0079
.022"x.028"	-10°	7°	7°D	R	41-2160-0079	40-2160-0079
				L	41-2260-0079	40-2260-0079

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



#### Type Andrews Convertible - HG .051\"/>

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.022\"x.028"	-10°	6°	0°	R	41-2160-0087	—
				L	41-2260-0087	—



#### Type Andrews Convertible - HG .045\"/>

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018\"x.025"	-10°	7°D	7°D	R	41-1160-0033	—
				L	41-1260-0033	—
.022\"x.028"	-10°	7°D	7°D	R	41-2160-0033	—
				L	41-2260-0033	—



#### Type Andrews Non-Convertible - HG .045\"/>

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018\"x.025"	-10°	6°D	0°	R	41-1160-0086	—
				L	41-1260-0086	—
.022\"x.028"	-10°	6°D	0°	R	41-2160-0086	—
				L	41-2260-0086	—



#### Edgewise Convertible - HG .045\"/>

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018\"x.025"	0°	0°	0°	R	41-1160-0077	—
				L	41-1260-0077	—
.022\"x.028"	0°	0°	0°	R	41-2160-0077	—
				L	41-2260-0077	—



#### Edgewise Non-Convertible - HG .045\"/>

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018\"x.025"	0°	0°	0°	R	41-1160-0084	—
				L	41-1260-0084	—
.022\"x.028"	0°	0°	0°	R	41-2160-0084	—
				L	41-2260-0084	—

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Maxillary Double Tubes 1<sup>st</sup> Molar



## Edgewise Convertible - HG .045"/Occlusal

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	0°	6°D	6°D	R	41-1160-0081	—
				L	41-1260-0081	—
.022"x.028"	0°	6°D	6°D	R	41-2160-0081	—
				L	41-2260-0081	—



## Edgewise Non-Convertible - HG .045"/Occlusal

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.022"x.028"	0°	6°D	0°	R	41-2160-0085	—
				L	41-2260-0085	—



## Type Planché Convertible - HG .045"/Gingival

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.022"x.028"	-14°	14°D	0°	R	41-2160-0070	—
				L	41-2260-0070	—

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Maxillary Triple Tubes 1<sup>st</sup> Molar



**MAXILLARY Triple Tubes 1st Molar**  
GET HG .045"/Occlusal



						P/N	
Dimensions	Torque	Rotation Slot	Rotation HG	R/L	Weldable	Bondable	
.022"x.028"	-20°	10°D	0°	R	41-2160-0500	—	
				L	41-2260-0500	—	
.022"x.028"	-20°	0°	0°	R	41-2160-0501	—	
				L	41-2260-0501	—	



**Type Roth Convertible - HG .045"/Occlusal**

						P/N	
Dimensions	Torque	Rotation Slot	Rotation HG	R/L	Weldable	Bondable	
.018"x.025"/ .018"x.025"	-14°	14°D	14°D	R	41-1160-0010	—	
				L	41-1260-0010	—	
.022"x.028"/ .018"x.025"	-14°	14°D	14°D	R	41-2160-0010	—	
				L	41-2260-0010	—	



**Type Roth Convertible - HG .045"/Occlusal**

						P/N	
Dimensions	Torque	Rotation Slot	Rotation HG	R/L	Weldable	Bondable	
.022"x.028"/ .022"x.028"	-14°	14°D	14°D	R	41-2160-0050	—	
				L	41-2260-0050	—	



**Type Roth Convertible - HG .045"/Occlusal**

						P/N	
Dimensions	Torque	Rotation Slot	Rotation HG	R/L	Weldable	Bondable	
.018"x.025"	-14°	0°	0°	R	41-1160-0051	—	
				L	41-1260-0051	—	
.022"x.028"	-14°	0°	0°	R	41-2160-0051	—	
				L	41-2260-0051	—	



**Type Roth Convertible - HG .051"/Occlusal**

						P/N	
Dimensions	Torque	Rotation Slot	Rotation HG	R/L	Weldable	Bondable	
.022"x.028"	-14°	14°D	14°D	R	41-2160-0052	—	
				L	41-2260-0052	—	

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Maxillary Triple Tubes 1<sup>st</sup> Molar



Type Roth Convertible - HG .051"/Occlusal

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.022"x.028"	-14°	0°	0°	R	41-2160-0053	—
				L	41-2260-0053	—



Type MBT\* Convertible F.E.O. .045"/Occlusal

Dimensions	Torque	Rotation Slot	Slot aux.	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-14°	10°	.018"x.025"	R	41-1160-0111	—
				L	41-1260-0111	—
.022"x.028"	-14°	10°	.018"x.025"	R	41-2160-0111	—
				L	41-2260-0111	—



Type Bio-Progressive Convertible - HG .045"/Gingival

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	0°	15°D	15°D	R	41-1160-0042	—
				L	41-1260-0042	—



Type Bio-Progressive Convertible - HG .045"/Gingival

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	0°	0°	0°	R	41-1160-0041	—
				L	41-1260-0041	—
.022"x.028"	0°	0°	0°	R	41-2160-0041	—
				L	41-2260-0041	—



Type Bio-Progressive Non-Convertible - HG .045"/Occlusal

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	0°	0°	0°	R	41-1160-0040	—
				L	41-1260-0040	—
.022"x.028"	0°	0°	0°	R	41-2160-0040	—
				L	41-2260-0040	—

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs


**Edgewise Convertible - HG .045"/Occlusal**

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	0°	6°D	0°	R	41-1160-0065	—
				L	41-1260-0065	—
.022"x.028"	0°	6°D	0°	R	41-2160-0065	—
				L	41-2260-0065	—


**Type Polar Convertible - HG .045"/Gingival**

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-20°	10°D	0°	R	41-1160-0027	—
				L	41-1260-0027	—


**Type Hilgers Convertible - HG .045"/Occlusal**

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-10°	15°D	0°	R	41-1160-0046	—
				L	41-1260-0046	—
.022"x.028"	-10°	15°D	0°	R	41-2160-0046	—
				L	41-2260-0046	—


**Type Melsen Convertible - HG .045"/Occlusal**

Dimensions	Torque	Rotation Slot	Rotation HG	R/L	P/N	
					Weldable	Bondable
.018"x.025"	-10°	8°D	0°	R	41-1160-0068	—
				L	41-1260-0068	—
.022"x.028"	-10°	8°D	0°	R	41-2160-0068	—
				L	41-2260-0068	—

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Mandibular Single Tubes 1<sup>st</sup> Molar



Convertible



## Type Roth

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-25°	0°	4°D	4.0	R	41-1460-0012	40-1460-0012
					L	41-1360-0012	40-1360-0012
.018"x.025"	-30°	0°	4°D	4.0	R	41-1460-0062	40-1460-0062
					L	41-1360-0062	40-1360-0062
.022"x.028"	-25°	0°	4°D	4.0	R	41-2460-0012	40-2460-0012
					L	41-2360-0012	40-2360-0012
.022"x.028"	-30°	0°	4°D	4.0	R	41-2460-0062	40-2460-0062
					L	41-2360-0062	40-2360-0062



## Type MBT\*

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-20°	0°	0°	3.8	R	41-1460-0022	40-1460-0119
					L	41-1360-0022	40-1360-0119
.018"x.025"	-10°	0°	0°	4.0	R	41-1460-0076	40-1460-0076
					L	41-1360-0076	40-1360-0076
.022"x.028"	-20°	0°	0°	3.8	R	41-2460-0119	40-2460-0119
					L	41-2360-0119	40-2360-0119
.022"x.028"	-10°	0°	0°	4.0	R	41-2460-0076	40-2460-0076
					L	41-2360-0076	40-2360-0076



## Type Alexander

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-22°	-6°	5°D	3.8	R	41-1460-0074	—
					L	41-1360-0074	—



## Type Andrews

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-25°	0°	0°	3.6	R	41-1460-0034	40-1460-0034
					L	41-1360-0034	40-1360-0034
.022"x.028"	-25°	0°	0°	3.6	R	41-2460-0034	40-2460-0034
					L	41-2360-0034	40-2360-0034

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

## Convertible



### Type Planché

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.022"x.028"	-25°	-6°	6°D	4.0	R	41-2460-0121	—
					L	41-2360-0121	—



### Edgewise

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.022"x.028"	0°	0°	6°D	4.0	R	—	40-2460-0092
					L	—	40-2360-0092



### Edgewise

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.018"x.025"	0°	0°	0°	4.0	R	-	40-1460-0093
					L	-	40-1360-0093
.022"x.028"	0°	0°	0°	4.0	R	41-2460-0093	40-2460-0093
					L	41-2360-0093	40-2360-0093



### Type Polar

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-22°	0°	4°D	4.0	R	—	40-2460-0029
					L	—	40-2360-0029

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Mandibular Single Tubes 2<sup>nd</sup> Molar



Non-Convertible




Type Roth

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-25°	0°	6°D	4.0	R	41-1470-0012	40-1470-0012
					L	41-1370-0012	40-1370-0012
.022"x.028"	-25°	0°	6°D	4.0	R	41-2470-0012	40-2470-0012
					L	41-2370-0012	40-2370-0012



Type Roth

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-30°	-1°	4°D	4.0	R	41-1470-0064	40-1470-0064
					L	41-1370-0064	40-1370-0064
.022"x.028"	-30°	-1°	4°D	4.0	R	41-2470-0064	40-2470-0064
					L	41-2370-0064	40-2370-0064



Type MBT\*

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-10°	0°	0°	3.2	R	41-1470-0022	40-1470-0022
					L	41-1370-0022	40-1370-0022
.022"x.028"	-10°	0°	0°	3.2	R	41-2470-0022	40-2470-0022
					L	41-2370-0022	40-2370-0022



Plan. Type

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.022"x.028"	-27°	-10°	6°D		R	41-2470-0012	—
					L	41-2370-0012	—



Alexander Type

Dimensions	Torque	Angulation	Rotation	Length mm	R/L	P/N	
						Weldable	Bondable
.022"x.028"	-27°	-10°	6°D		R	—	40-2470-0074
					L	—	40-2370-0074

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

## Non-Convertible



### Type Andrews

							P/N	
Dimensions	Torque	Angulation	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-25°	0°	0°	4.0	R	—	40-1470-0034	
					L	—	40-1370-0034	
.022"x.028"	-25°	0°	0°	4.0	R	—	40-2470-0034	
					L	—	40-2370-0034	
.018"x.025"	-25°	0°	4°D	3.0	R	—	40-1470-0049	
					L	—	40-1370-0049	
.022"x.028"	25°	0°	4°D	3.0	R	—	40-2470-0049	
					L	—	40-2370-0049	



### Type Andrews Wide Entry (Peucat)

							P/N	
Dimensions	Torque	Angulation	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-25°	0°	4°D	3.5	R	41-1470-0018	40-1470-0018	
					L	41-1370-0018	40-1370-0018	
.022"x.028"	-25°	0°	4°D	3.5	R	41-2470-0018	40-2470-0018	
					L	41-2370-0018	40-2370-0018	



### Type Polar

							P/N	
Dimensions	Torque	Angulation	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-25°	5°	7°D	4	R	—	40-1470-0029	
					L	—	40-1370-0029	



### Edgewise

							P/N	
Dimensions	Torque	Angulation	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	0°	0°	0°	4.6	R	—	40-1470-0094	
					L	—	40-1370-0094	
.022"x.028"	0°	0°	0°	4.6	R	41-2470-0094	40-2470-0094	
					L	41-2370-0094	40-2370-0094	



### Edgewise

							P/N	
Dimensions	Torque	Angulation	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	0°	0°	0°	4.0	R	41-1470-0095	40-1470-0095	
					L	41-1370-0095	40-1370-0095	
.022"x.028"	0°	0°	0°	4.0	R	41-2470-0095	40-2470-0095	
					L	41-2370-0095	40-2370-0095	

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Mandibular Single Tubes 2<sup>nd</sup> Molar



Non-Convertible



Edgewise

							P/N	
Dimensions	Torque	Angulation	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	0°	0°	6°D	4.0	R	—	40-1470-0096	
					L	—	40-1370-0096	
.022"x.028"	0°	0°	6°D	4.0	R	41-2470-0096	40-2470-0096	
					L	41-2370-0096	40-2370-0096	



Edgewise

							P/N	
Dimensions	Torque	Angulation	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	0°	0°	6°D	3.2	R	—	40-1470-0097	
					L	—	40-1370-0097	
.022"x.028"	0°	0°	6°D	3.2	R	41-2470-0096	40-2470-0097	
					L	41-2370-0096	40-2370-0097	



Edgewise

							P/N	
Dimensions	Torque	Angulation	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	0°	0°	0°	3.2	R	—	40-1470-0098	
					L	—	40-1370-0098	
.022"x.028"	0°	0°	0°	3.2	R	41-2470-0098	40-2470-0098	
					L	41-2370-0098	40-2370-0098	



Small

							P/N	
Dimensions	Torque	Angulation	Rotation	Length mm	R/L	Weldable	Bondable	
.018"x.025"	-25°	0°	4°D	3.0	R	—	40-1470-0048	
					L	—	40-1370-0048	
.022"x.028"	-25°	0°	4°D	3.0	R	—	40-2470-0048	
					L	—	40-2370-0048	



Mini Eruption Tube for 2<sup>nd</sup> Molar

							P/N	
Dimensions	Torque	Angulation	Rotation	Length mm	R/L	Without hook	With hook	
.018"x.025"	-25°	0°	0°	2.0	R	40-1470-0112	-	
					L	40-1370-0112	-	
.022"x.028"	-25°	0°	0°	2.0	R	40-2470-0112	40-2470-0019	
					L	40-2370-0112	40-2370-0019	

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Mandibular Double Tubes 1<sup>st</sup> Molar



Convertible with Tube Bumper



Type BiDimensional .051"/Gingival

Dimensions	Torque	Angulation Tube	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.022"x.028"	0°	0°	0°	0°	R	41-2460-0067	—
					L	41-2360-0067	—



Type Cetlin .045"/Gingival

Dimensions	Torque	Angulation Tube	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-25°	0°	4°D	5M	R	41-1460-0037	—
					L	41-1360-0037	—
.022"x.028"	-25°	0°	4°D	5M	R	41-2460-0037	—
					L	41-2360-0037	—

Convertible



Type Roth Convertible

Dimensions	Torque	Angulation Tube	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.018"x.025" -	25°	0°	4°D	0°	R	41-1460-0011	40-1460-0011
					L	41-1360-0011	40-1360-0011
.022"x.028"	-25°	0°	4°D	0°	R	41-2460-0011	40-2460-0011
					L	41-2360-0011	40-2360-001



Type Roth Occlusal Hook

Dimensions	Torque	Angulation Tube	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.022"x.028"	-30°	-1°	4°D	0°	R	41-2460-0058	40-2460-0058
					L	41-2360-0058	40-2360-0058

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



# Mandibular Double Tubes 1<sup>st</sup> Molar



Convertible



## Type Roth Gingival Hook

Dimensions	Torque	Angulation Tube	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.022"x.028"	-30°	-1°	4°D	0°	R	41-2460-0060	40-2460-0060
					L	41-2360-0060	40-2360-0060




## Extra Torque

Dimensions	Torque	Angulation Tube	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-10°	0°	4°D	0°	R	41-1460-0070	—
					L	41-1360-0070	—
.022"x.028"	-10°	0°	4°D	0°	R	41-2460-0070	—
					L	41-2360-0070	—



## Extra Torque Wide Base

Dimensions	Torque Slot & Aux.	Angulation Tube	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-10°	0°	4°D	0°	R	—	40-1460-0031
					L	—	40-1360-0031



## Type Bio-Progressive

Dimensions	Torque Slot	Torque Aux.	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.018"x.025"	0°	0°	0°	0°	UNI	41-1660-0045	40-1660-0045
.022"x.028"	0°	0°	0°	0°	UNI	41-2660-0045	40-2660-0045



## Type Polar

Dimensions	Torque Slot	Torque Aux.	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-25°	0°	4°D	0°	R	41-1460-0028	—
					L	41-1360-0028	—

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

### Convertible



#### Type Hilgers

Dimensions	Torque Slot	Torque Aux.	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-27°	0°	12°D	0°	R	41-1460-0047	—
					L	41-1360-0047	—
.022"x.028"	-27°	0°	12°D	0°	R	41-2460-0047	—
					L	41-2360-0047	—



#### Type Melsen

Dimensions	Torque Slot	Torque Aux.	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.018"x.025"	-27°	0°	8°D	0°	R	41-1460-0069	—
					L	41-1360-0069	—
.022"x.028"	-27°	0°	8°D	0°	R	41-2460-0069	—
					L	41-2360-0069	—

### Non- Convertible



#### Type Roth

Dimensions	Torque	Angulation Tube	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.022"x.028"	-30°	-1°	4°D	0°	R	41-2460-0059	—
					L	41-2360-0059	—








#### Type Bio-Progressive





Dimensions	Torque Slot	Torque Aux.	Rotation Slot	Rotation Aux.	R/L	P/N	
						Weldable	Bondable
.018"x.025"	0°	0°	0°	0°	UNI	41-1660-0044	40-1660-0044
.022"x.028"	0°	0°	0°	0°	UNI	41-2660-0044	40-2660-0044

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

# WELDABLE Lingual Attachments

Lingual Sheats		P/N
	Standard	43-0250-0000
	Cetlin left 8° offset, V-Stop, Hook	43-027L-0000
	Cetlin right 8° offset, V-Stop, Hook	43-027R-0000
	With left Hook, V-Stop	43-028L-0000
	With right Hook, V-Stop	43-028R-0000
	Without V stop	43-0260-0000
	With rotation wings and stops, Upper left and lower right	43-042L-0000
	With rotation wings and stops, Upper right and lower left	43-042R-0000
	Curved Without V-Stops	43-0430-0000

Seating Lugs		P/N
	Flat	43-0340-0000
	Curved	43-0350-0000

Cleat Lingual		P/N
	Small	43-0210-0000
	Medium	43-0220-0000
	Large	43-0230-0000
	Large Optional	43-0240-0000




# VIPER Tubes

## Buccal Single Non-Convertible

### VIPER™ LOW PROFILE BUCCAL TUBES

Innovative technology and sleek, stylish comfort

VIPER™ Low Profile Buccal Tubes Sleek and stylish, G&H®'s Viper™ buccal tube incorporates the latest design technology with a growing industry demand for comfortable, easy-to-use appliances that simplify treatment and reduce valuable chair time.



#### MAXILLARY 1<sup>st</sup> & 2<sup>nd</sup> Molar

Torque	Offset	Rx	M/D	R/L	ID	Bondable			
						.018	.022	.018	.022
0°	0°	STD	4.2 mm	R	◆	NS3B-821E	NS3B-221E	NS3B-821I	NS3B-221I
				L	◆	NS3B-822E	NS3B-222E	NS3B-822I	NS3B-222I
-10°	8°	Roth	4.2 mm	R	●	NS3B-825E	NS3B-225E	NS3B-825I	NS3B-225I
				L	●	NS3B-826E	NS3B-226E	NS3B-826I	NS3B-226I
-14°	10°	MBT	4.2 mm	R	●	NS3B-833E	NS3B-233E	NS3B-833I	NS3B-233I
				L	●	NS3B-834E	NS3B-234E	NS3B-834I	NS3B-234I

#### MANDIBULAR 1<sup>st</sup> & 2<sup>nd</sup> Molar

Torque	Offset	Rx	M/D	R/L	ID	Bondable			
						.018	.022	.018	.022
0°	0°	STD	4.2 mm	R	◆	NS3B-822E	NS3B-222E	NS3B-822I	NS3B-222I
				L	◆	NS3B-821E	NS3B-221E	NS3B-821I	NS3B-221I
-25°	5°	Roth	4.2 mm	R	■	NS3B-827E	NS3B-227E	NS3B-827I	NS3B-227I
				L	■	NS3B-828E	NS3B-228E	NS3B-828I	NS3B-228I
-20°	0°	MBT	4.2 mm	R	▲	NS3B-809E	NS3B-209E	NS3B-809I	NS3B-209I
				L	▲	NS3B-810E	NS3B-210E	NS3B-810I	NS3B-210I

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

# Evolve Flow Bondable Tubes

## Positioning Guide


Occlusal/gingival positioning guides ensure a secure grip and ease of placement

## Distal Notch

Evolve Flow features a distal notch for ligation and tieback of auxiliaries

## Contoured Base

Anatomically contoured base with occlusal indent ensures accurate fit, whilst the 80 gauge mesh pad ensures superior bond strength



## Low profile

The contoured surfaces are low profile and tapered for patient comfort

## Low profile hook

Low profile hooks guarantee patient comfort and allow for easy placements of attachments

## The Funnelled Entrance

The Evolve Flow funnelled entrance allows for easy wire insertion

Roth .022	1 <sup>st</sup> Molars	2 <sup>nd</sup> Molars
Upper Left	DB07-1350UL	DB07-1371UL
Upper Right	DB07-1350UR	DB07-1371UR
Lower Left	DB07-1350LL	DB07-1371LL
Lower Right	DB07-1350LR	DB07-1371LR

MBT .022	1 <sup>st</sup> Molars	2 <sup>nd</sup> Molars
Upper Left	DB07-1360UL	DB07-1372UL
Upper Right	DB07-1360UR	DB07-1372UR
Lower Left	DB07-1360LL	DB07-1372LL
Lower Right	DB07-1360LR	DB07-1372LR

Roth .018	1 <sup>st</sup> Molars	2 <sup>nd</sup> Molars
Upper Left	DB07-1370UL	DB07-1373UL
Upper Right	DB07-1370UR	DB07-1373UR
Lower Left	DB07-1370LL	DB07-1373LL
Lower Right	DB07-1370LR	DB07-1373LR

MBT .018	1 <sup>st</sup> Molars	2 <sup>nd</sup> Molars
Upper Left	DB07-1362UL	DB07-1374UL
Upper Right	DB07-1362UR	DB07-1374UR
Lower Left	DB07-1362LL	DB07-1374LL
Lower Right	DB07-1362LR	DB07-1374LR

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

