



**Kleef- & Weldable tubes/
Tubes à coller & à souder**
TUBES



SELF-LIGATING Tubes



Experience Self-ligating tubes

MAXILLARY

| Tooth 6 & 7 | Torque | Angulation | Rotation | .018" | | .022" | |
|---|--------|------------|----------|--------------|--------------|--------------|--------------|
| | | | | R | L | R | L |
| Experience SL tube, Roth | -10° | 0° | 14°D | 20-1160-0010 | 20-1260-0010 | 20-2160-0010 | 20-2260-0010 |
| Experience SL tube, Small Base Roth | -10° | 0° | 14°D | - | - | 20-2160-S010 | 20-2260-S010 |
| Experience SL tube, MBT | -14° | 0° | 10°D | 20-1160-0020 | 20-1260-0020 | 20-2160-0020 | 20-2260-0020 |
| Experience SL tube, Small Base MBT | -14° | 0° | 10°D | 20-1160-S020 | 20-1260-S020 | 20-2160-S020 | 20-2260-S020 |
| Experience double SL tube, Extra Torque | -20° | 0° | 20°D | 20-1160-0030 | 20-1260-0030 | 20-2160-0030 | 20-2260-0030 |

MANDIBULAR

| Tooth 6 & 7 | Torque | Angulation | Rotation | .018" | | .022" | |
|-------------------------------------|--------|------------|----------|--------------|--------------|--------------|--------------|
| | | | | R | L | R | L |
| Experience SL tube, Roth | -25° | 0° | 4°D | 20-1460-0010 | 20-1360-0010 | 20-2460-0010 | 20-2360-0010 |
| Experience SL tube, Small Base Roth | -25° | 0° | 4°D | - | - | 20-2460-S010 | 20-2360-S010 |
| Experience SL tube, MBT | -20° | 0° | 0°D | 20-1460-0020 | 20-1360-0020 | 20-2460-0020 | 20-2360-0020 |
| Experience SL tube, Small Base MBT | -20° | 0° | 0°D | 20-1460-S020 | 20-1360-S020 | 20-2460-S020 | 20-2360-S020 |
| Experience SL tube, Extra Torque | -14° | 0° | 4°D | - | - | 20-2460-0032 | 20-2360-0032 |

MANDIBULAR

| Tooth 6 & 7 | Torque | Angulation | Rotation | .018" | | .022" | |
|---|--------|------------|----------|--------------|--------------|--------------|--------------|
| | | | | R | L | R | L |
| Experience DOUBLE SL tube, Extra Torque | -10° | 0° | 4°D | 20-1460-0030 | 20-1360-0030 | 20-2460-0030 | 20-2360-0030 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



LEGEND LP Tubes

Conventional tubes

| Tooth | Torque | Angulation | Rotation | .018" | | .022" | |
|-----------|--------|------------|----------|--------------|--------------|--------------|--------------|
| | | | | R | L | R | L |
| U6 , Roth | -14° | 0° | 14° | 40-1160-0114 | 40-1260-0114 | 40-2160-0114 | 40-2260-0114 |
| U7, Roth | -14° | 0° | 14° | 40-1170-0114 | 40-1270-0114 | 40-2170-0114 | 40-2270-0114 |
| L6, Roth | -25° | 0° | 4° | 40-1460-0114 | 40-1360-0114 | 40-2460-0114 | 40-2360-0114 |
| L7, Roth | -25° | 0° | 6° | 40-1470-0114 | 40-1370-0114 | 40-2470-0114 | 40-2370-0114 |
| U6 , MBT | -14° | 0° | 10° | 40-1160-0106 | 40-1260-0106 | 40-2160-0106 | 40-2260-0106 |
| U6 , MBT | -14° | 0° | 0° | 40-1160-0107 | 40-1260-0107 | 40-2160-0107 | 40-2260-0107 |
| U7, MBT | -14° | 0° | 10° | 40-1170-0106 | 40-1270-0106 | 40-2170-0106 | 40-2270-0106 |
| U7, MBT | -14° | 0° | 0° | 40-1170-0107 | 40-1270-0107 | 40-2170-0107 | 40-2270-0107 |
| L6, MBT | -20° | 0° | 0° | 40-1460-0106 | 40-1360-0106 | 40-2460-0106 | 40-2360-0106 |
| L7, MBT | -10° | 0° | 0° | 40-1470-0106 | 40-1370-0106 | 40-2470-0106 | 40-2370-0106 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

Kleef- & Weldable TUBES à coller & à souder



LEGEND LE Tubes



LEGEND LE Tubes (Large Entrance)

Second Molar Bondable

Non-convertible single tubes for second molars, the funnel entry of these tubes facilitates access to second molars and its anatomical base fits the teeth perfectly making positioning easy.

MAXILLARY

Type Roth

| Dimensions | Torque | Rotation | P/N | |
|-------------|--------|----------|--------------|--------------|
| | | | R | L |
| .018"x.025" | -14° | 14°D | 40-1170-0099 | 40-1270-0099 |
| .022"x.028" | -14° | 14°D | 40-2170-0099 | 40-2270-0099 |

Type Andrews

| Dimensions | Torque | Rotation | P/N | |
|-------------|--------|----------|--------------|--------------|
| | | | R | L |
| .018"x.025" | -9° | 7°D | 40-1170-0100 | 40-1270-0100 |
| .022"x.028" | -9° | 7°D | 40-2170-0100 | 40-2270-0100 |

MANDIBULAR

Type Roth

| Dimensions | Torque | Rotation | P/N | |
|-------------|--------|----------|--------------|--------------|
| | | | R | L |
| .018"x.025" | -30° | 6°D | 40-1470-0099 | 40-1370-0099 |
| .022"x.028" | -30° | 6°D | 40-2470-0099 | 40-2370-0099 |


Type Andrews

| Dimensions | Torque | Rotation | P/N | |
|-------------|--------|----------|--------------|--------------|
| | | | R | L |
| .018"x.025" | -25° | 0° | 40-1470-0100 | 40-1370-0100 |
| .022"x.028" | -25° | 0° | 40-2470-0100 | 40-2370-0100 |

Type MBT*

| Dimensions | Torque | Rotation | P/N | |
|-------------|--------|----------|--------------|--------------|
| | | | R | L |
| .018"x.025" | -10° | 0° | 40-1470-0101 | 40-1370-0101 |
| .022"x.028" | -10° | 0° | 40-2470-0101 | 40-2370-0101 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



AXCESS Tubes



AXCESS Tubes

These entry point buccal tubes are manufactured in one-piece using MIM (Metal Injection Molded) technology. AXCESS tubes are available in bondable and weldable versions.

MAXILLARY

Double Non-Convertible HG .045" 1st Molar

| Dimensions | Torque | Rotation | R/L | P/N | |
|-------------|--------|----------|-----|--------------|----------|
| | | | | Weldable | Bondable |
| .018"x.025" | -10° | 6°D | R | 41-1160-0023 | — |
| | | | L | 41-1260-0023 | — |
| .022"x.028" | -10° | 6°D | R | 41-2160-0023 | — |
| | | | L | 41-2260-0023 | — |
| .018"x.025" | 0° | 6°D | R | 41-1160-0024 | — |
| | | | L | 41-1260-0024 | — |
| .022"x.028" | 0° | 6°D | R | 41-2160-0024 | — |
| | | | L | 41-2260-0024 | — |

Single Non-Convertible 1st Molar

| Dimensions | Torque | Rotation | R/L | P/N | |
|-------------|--------|----------|-----|--------------|--------------|
| | | | | Weldable | Bondable |
| .022"x.028" | -10° | 6°D | R | 41-2160-0025 | 40-2160-0025 |
| | | | L | 41-2260-0025 | 40-2260-0025 |
| .018"x.025" | 0° | 6°D | R | 41-1160-0026 | 40-1160-0026 |
| | | | L | 41-1260-0026 | 40-1260-0026 |
| .022"x.028" | 0° | 6°D | R | — | 40-2160-0026 |
| | | | L | — | 40-2260-0026 |

Single Non-Convertible 2nd Molar

| Dimensions | Torque | Rotation | R/L | P/N | |
|-------------|--------|----------|-----|----------|--------------|
| | | | | Weldable | Bondable |
| .018"x.025" | -10° | 6°D | R | — | 40-1170-0088 |
| | | | L | — | 40-1270-0088 |
| .018"x.025" | 0° | 0° | R | — | 40-1170-0090 |
| | | | L | — | 40-1270-0090 |
| .018"x.025" | -10° | 0° | R | — | 40-1170-0091 |
| | | | L | — | 40-1270-0091 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



AXCESS Tubes



AXCESS Tubes

Single Non-Convertible 2nd Molar

| Dimensions | Torque | Rotation | R/L | P/N | |
|-------------|--------|----------|-----|--------------|--------------|
| | | | | Weldable | Bondable |
| .022"x.028" | -10° | 6°D | R | 41-2170-0025 | 40-2170-0025 |
| | | | L | 41-2270-0025 | 40-2270-0025 |
| .022"x.028" | -10° | 6°D | R | — | 40-2170-0088 |
| | | | L | — | 40-2270-0088 |
| .022"x.028" | 0° | 0° | R | — | 40-2170-0090 |
| | | | L | — | 40-2270-0090 |
| .022"x.028" | -10° | 0° | R | — | 40-2170-0091 |
| | | | L | — | 40-2270-0091 |

MANDIBULAR

Single Non-Convertible 1st Molar

| Dimensions | Torque | Rotation | R/L | P/N | |
|-------------|--------|----------|-----|--------------|--------------|
| | | | | Weldable | Bondable |
| .022"x.028" | -25° | 6°D | R | 41-2460-0025 | 40-2460-0025 |
| | | | L | 41-2360-0025 | 40-2360-0025 |
| .018"x.025" | -25° | 6°D | R | 41-1460-0025 | - |
| | | | L | 41-1360-0025 | - |
| .022"x.028" | 0° | 6°D | R | 41-2460-0026 | 40-2460-0026 |
| | | | L | 41-2360-0026 | 40-2360-0026 |
| .018"x.025" | 0° | 6°D | R | 41-1460-0026 | 40-1460-0026 |
| | | | L | 41-1360-0026 | 40-1360-0026 |

Single Non-Convertible 2nd Molar

| Dimensions | Torque | Rotation | R/L | P/N | |
|-------------|--------|----------|-----|--------------|--------------|
| | | | | Weldable | Bondable |
| .018"x.025" | 0° | 0° | R | — | 40-1470-0090 |
| | | | L | — | 40-1370-0090 |
| .018"x.025" | -25° | 0° | R | — | 40-1470-0091 |
| | | | L | — | 40-1370-0091 |
| .022"x.028" | -25° | 6°D | R | 41-2470-0025 | 40-2470-0025 |
| | | | L | 41-2370-0025 | 40-2370-0025 |
| .022"x.028" | -25° | 6°D | R | — | 40-2470-0089 |
| | | | L | — | 40-2370-0089 |
| .022"x.028" | 0° | 0° | R | — | 40-2470-0090 |
| | | | L | — | 40-2370-0090 |
| .022"x.028" | -25° | 6°D | R | — | 40-2470-0091 |
| | | | L | — | 40-2370-0091 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Maxillary Single Tubes Convertible



Type Roth, 1st Molar

| | | | | | | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -14° | 14°D | 4.0 | R | 41-1160-0061 | 40-1160-0061 | |
| | | | | L | 41-1260-0061 | 40-1260-0061 | |
| .022"x.028" | -14° | 14°D | 4.0 | R | 41-2160-0061 | 40-2160-0061 | |
| | | | | L | 41-2260-0061 | 40-2260-0061 | |



Type Roth, 1st Molar

| | | | | | | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -10° | 14°D | 4.0 | R | 41-1160-0012 | 40-1160-0012 | |
| | | | | L | 41-1260-0012 | 40-1260-0012 | |
| .022"x.028" | -10° | 14°D | 4.0 | R | 41-2160-0012 | 40-2160-0012 | |
| | | | | L | 41-2260-0012 | 40-2260-0012 | |



Type MBT*, 1st Molar

| | | | | | | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -14° | 10°D | 3.8 | R | 41-1160-0022 | 40-1160-0022 | |
| | | | | L | 41-1260-0022 | 40-1260-0022 | |
| .022"x.028" | -14° | 10°D | 3.8 | R | 41-2160-0022 | 40-2160-0022 | |
| | | | | L | 41-2260-0022 | 40-2260-0022 | |



Type MBT*, 1st Molar

| | | | | | | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -14° | 10°D | 3.8 | R | 41-1160-0110 | 40-1160-0110 | |
| | | | | L | 41-1260-0110 | 40-1260-0110 | |
| .022"x.028" | -14° | 10°D | 3.8 | R | 41-2160-0110 | 40-2160-0110 | |
| | | | | L | 41-2260-0110 | 40-2260-0110 | |



Type Andrews, 1st Molar

| | | | | | | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -10° | 7°D | 3.6 | R | 41-1160-0034 | 40-1160-0034 | |
| | | | | L | 41-1260-0034 | 40-1260-0034 | |
| .022"x.028" | -10° | 7°D | 3.6 | R | 41-2160-0034 | 40-2160-0034 | |
| | | | | L | 41-2260-0034 | 40-2260-0034 | |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Maxillary Single Tubes Non-Convertible



GET 1st Molar



| Dimensions | Torque | Rotation | Length mm | R/L | P/N | |
|-------------|--------|----------|-----------|-----|----------|--------------|
| | | | | | Weldable | Bondable |
| .022"x.028" | -20° | 10°D | 4.0 | R | — | 40-2160-0510 |
| | | | | L | — | 40-2260-0510 |
| .022"x.028" | -20° | 0° | 4.0 | R | — | 40-2160-0511 |
| | | | | L | — | 40-2260-0511 |



GET 2nd Molar



| Dimensions | Torque | Rotation | Length mm | R/L | P/N | |
|-------------|--------|----------|-----------|-----|----------|--------------|
| | | | | | Weldable | Bondable |
| .022"x.028" | -20° | 0° | 4.0 | R | — | 40-2170-0120 |
| | | | | L | — | 40-2270-0120 |



Type Roth Distal Extension, 2nd Molar

| Dimensions | Torque | Rotation | Length mm | R/L | P/N | |
|-------------|--------|----------|-----------|-----|----------|--------------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -14° | 14°D | 4.0 | R | — | 40-1170-0013 |
| | | | | L | — | 40-1270-0013 |
| .022"x.028" | -14° | 14°D | 4.0 | R | — | 40-2170-0013 |
| | | | | L | — | 40-2270-0013 |



Type Roth Distal Extension, 2nd Molar

| Dimensions | Torque | Rotation | Length mm | R/L | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|
| | | | | | Weldable | Bondable |
| .022"x.028" | -10° | 14°D | 3.0 | R | 41-2170-0039 | 40-2170-0039 |
| | | | | L | 41-2270-0039 | 40-2270-0039 |



Type Roth, 2nd Molar

| Dimensions | Torque | Rotation | Length mm | R/L | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -14° | 14°D | 4.2 | R | — | 40-1170-0014 |
| | | | | L | — | 40-1270-0014 |
| .022"x.028" | -14° | 14°D | 4.2 | R | 41-2170-0014 | 40-2170-0014 |
| | | | | L | 41-2270-0014 | 40-2270-0014 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs




Maxillary Single Tubes Non-Convertible



Type Roth, 2nd Molar

| Dimensions | Torque | Rotation | Length mm | R/L | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -14° | 0° | 4.0 | R | — | 40-1170-0063 |
| | | | | L | — | 40-1270-0063 |
| .022"x.028" | -14° | 0° | 4.0 | R | 41-2170-0063 | 40-2170-0063 |
| | | | | L | 41-2270-0063 | 40-2270-0063 |



Type Roth, 2nd Molar

| Dimensions | Torque | Rotation | Length mm | R/L | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -14° | 14°D | 4.0 | R | 41-1170-0012 | 40-1170-0012 |
| | | | | L | 41-1270-0012 | 40-1270-0012 |
| .022"x.028" | -14° | 14°D | 4.0 | R | 41-2170-0012 | 40-2170-0012 |
| | | | | L | 41-2270-0012 | 40-2270-0012 |



Type MBT*, 2nd Molar

| Dimensions | Torque | Rotation | Length mm | R/L | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -14° | 10°D | 4.2 | R | 41-1170-0022 | 40-1170-0022 |
| | | | | L | 41-1270-0022 | 40-1270-0022 |
| .022"x.028" | -14° | 10°D | 4.2 | R | 41-2170-0022 | 40-2170-0022 |
| | | | | L | 41-2270-0022 | 40-2270-0022 |



Type Alexander Distal Extension, 2nd Molar

| Dimensions | Torque | Rotation | Length mm | R/L | P/N | |
|-------------|--------|----------|-----------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .022"x.028" | -10° | 7°D | 4.0 | R | 41-2170-0075 | — |
| | | | | L | 41-2270-0075 | — |



Type Andrews, 2nd Molar


| Dimensions | Torque | Rotation | Length mm | R/L | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -10° | 7°D | 4.0 | R | 41-1170-0034 | 40-1170-0034 |
| | | | | L | 41-1270-0034 | 40-1270-0034 |
| .022"x.028" | -10° | 7°D | 4.0 | R | 41-2170-0034 | 40-2170-0034 |
| | | | | L | 41-2270-0034 | 40-2270-0034 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Type Andrews, 2nd Molar

| | | | | | | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -10° | 7°D | 3.0 | R | 41-1170-0048 | 40-1170-0048 | |
| | | | | L | 41-1270-0048 | 40-1270-0048 | |
| .022"x.028" | -10° | 7°D | 3.0 | R | 41-2170-0048 | 40-2170-0048 | |
| | | | | L | 41-2270-0048 | 40-2270-0048 | |



Type Andrews Wide Base, 2nd Molar

| | | | | | | P/N | |
|-------------|--------|----------|-----------|-----|----------|--------------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -10° | 7°D | 3.0 | R | — | 40-1170-0049 | |
| | | | | L | — | 40-1270-0049 | |
| .022"x.028" | -10° | 7°D | 3.0 | R | — | 40-2170-0049 | |
| | | | | L | — | 40-2270-0049 | |




Type Andrews Wide Entry, 2nd Molar

| | | | | | | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -10° | 7°D | 3.5 | R | 41-1170-0018 | 40-1170-0018 | |
| | | | | L | 41-1270-0018 | 40-1270-0018 | |
| .022"x.028" | -10° | 7°D | 3.5 | R | 41-2170-0018 | 40-2170-0018 | |
| | | | | L | 41-2270-0018 | 40-2270-0018 | |



Type Planché, 2nd Molar

| | | | | | | P/N | |
|---------------|--------|----------|-----------|-----|--------------|----------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Weldable | Bondable | |
| .022"x.028" - | 20° | 14°D | 4.2 | R | 41-2170-0071 | — | |
| | | | | L | 41-2270-0071 | — | |



Type Polar, 2nd Molar

| | | | | | | P/N | |
|-------------|--------|----------|-----------|-----|----------|--------------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -25° | 6°D | 4.0 | R | — | 40-1170-0029 | |
| | | | | L | — | 40-1270-0029 | |



Mini Eruption Tube, 2nd Molar

| | | | | | | P/N | |
|-------------|--------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Rotation | Length mm | R/L | Without hook | With hook | |
| .018"x.025" | -10° | 0°D | 2.0 | R | 40-1170-0112 | 40-1170-0019 | |
| | | | | L | 40-1270-0112 | 40-1270-0019 | |
| .022"x.028" | -10° | 0°D | 2.0 | R | 40-2170-0112 | 40-2170-0019 | |
| | | | | L | 40-2270-0112 | 40-2270-0019 | |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Maxillary Double Tubes 1st Molar Convertible



Type Roth Convertible No HG

| Dimensions | Torque | Rotation Slot & Aux. | R/L | P/N | |
|---------------------------|--------|----------------------|-----|--------------|--------------|
| | | | | Weldable | Bondable |
| .018"x.025" / .018"x.025" | -14° | 14°D | R | 41-1160-0011 | 40-1160-0011 |
| | | | L | 41-1260-0011 | 40-1260-0011 |
| .022"x.028" / .018"x.025" | -14° | 14°D | R | 41-2160-0011 | 40-2160-0011 |
| | | | L | 41-2260-0011 | 40-2260-0011 |



Type Roth Convertible No HG

| Dimensions | Torque | Rotation Slot | R/L | P/N | |
|---------------------------|--------|---------------|-----|--------------|--------------|
| | | | | Weldable | Bondable |
| .018"x.025" / .018"x.025" | -14° | 14°D | R | 41-1160-0055 | 40-1160-0055 |
| | | | L | 41-1260-0055 | 40-1260-0055 |
| .022"x.028" / .022"x.028" | -14° | 14°D | R | 41-2160-0055 | 40-2160-0055 |
| | | | L | 41-2260-0055 | 40-2260-0055 |



Type MBT* Convertible No HG


| Dimensions | Torque | Rotation Slot | R/L | P/N | |
|-------------|--------|---------------|-----|--------------|--------------|
| | | | | Weldable | Bondable |
| .018"x.025" | -14° | 10°D | R | 41-1160-0104 | 40-1160-0104 |
| | | | L | 41-1260-0104 | 40-1260-0104 |
| .022"x.028" | -14° | 10°D | R | 41-2160-0104 | 40-2160-0104 |
| | | | L | 41-2260-0104 | 40-2260-0104 |



Type MBT* Convertible No HG


| Dimensions | Torque | Rotation Slot | R/L | P/N | |
|-------------|--------|---------------|-----|--------------|--------------|
| | | | | Weldable | Bondable |
| .018"x.025" | -14° | 0°D | R | 41-1160-0105 | 40-1160-0105 |
| | | | L | 41-1260-0105 | 40-1260-0105 |
| .022"x.028" | -14° | 0°D | R | 41-2160-0105 | 40-2160-0105 |
| | | | L | 41-2260-0105 | 40-2260-0105 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Extra Torque Convertible No HG

| Dimensions | Torque | Rotation Slot | R/L | P/N | |
|-------------|--------|---------------|-----|--------------|----------|
| | | | | Weldable | Bondable |
| .022"x.028" | -20° | 20°D | R | 41-2160-0030 | - |
| | | | L | 41-2260-0030 | - |



Type Roth Convertible - HG .045"/Gingival with HG

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -14° | 14°D | 14°D | R | 41-1160-0054 | — |
| | | | | L | 41-1260-0054 | — |
| .022"x.028" | -14° | 14°D | 14°D | R | 41-2160-0054 | — |
| | | | | L | 41-2260-0054 | — |



Type MBT* Convertible - HG .045"/Occlusal With HG

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -14° | 10°D | 0° | R | 41-1160-0021 | — |
| | | | | L | 41-1260-0021 | — |
| .022"x.028" | -14° | 10°D | 0° | R | 41-2160-0021 | — |
| | | | | L | 41-2260-0021 | — |



Type MBT* Convertible - HG .045"/Gingival With HG

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .022"x.028" | -14° | 10°D | 10°D | R | 41-2160-0020 | - |
| | | | | L | 41-2260-0020 | - |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Maxillary Double Tubes 1st Molar



Type BiDimensional HG Convertible .051"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .022"x.028" | 0° | 6°D | 0° | R | 41-2160-0080 | — |
| | | | | L | 41-2260-0080 | — |



Type BiDimensional HG Convertible .051"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .022"x.028" | 0° | 14°D | 7°D | R | 41-2160-0083 | — |
| | | | | L | 41-2260-0083 | — |



Type Alexander Convertible - HG .045"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -10° | 14°D | 7°D | R | 41-1160-0073 | — |
| | | | | L | 41-1260-0073 | — |
| .022"x.028" | -10° | 14°D | 7°D | R | 41-2160-0073 | — |
| | | | | L | 41-2260-0073 | — |



Type Andrews Convertible - HG .045"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|--------------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -10° | 7° | 7°D | R | 41-1160-0079 | 40-1160-0079 |
| | | | | L | 41-1260-0079 | 40-1260-0079 |
| .022"x.028" | -10° | 7° | 7°D | R | 41-2160-0079 | 40-2160-0079 |
| | | | | L | 41-2260-0079 | 40-2260-0079 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs




Type Andrews Convertible - HG .051\"/>

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|--------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .022\"x.028" | -10° | 6° | 0° | R | 41-2160-0087 | — |
| | | | | L | 41-2260-0087 | — |



Type Andrews Convertible - HG .045\"/>

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|--------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018\"x.025" | -10° | 7°D | 7°D | R | 41-1160-0033 | — |
| | | | | L | 41-1260-0033 | — |
| .022\"x.028" | -10° | 7°D | 7°D | R | 41-2160-0033 | — |
| | | | | L | 41-2260-0033 | — |



Type Andrews Non-Convertible - HG .045\"/>

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|--------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018\"x.025" | -10° | 6°D | 0° | R | 41-1160-0086 | — |
| | | | | L | 41-1260-0086 | — |
| .022\"x.028" | -10° | 6°D | 0° | R | 41-2160-0086 | — |
| | | | | L | 41-2260-0086 | — |



Edgewise Convertible - HG .045\"/>

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|--------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018\"x.025" | 0° | 0° | 0° | R | 41-1160-0077 | — |
| | | | | L | 41-1260-0077 | — |
| .022\"x.028" | 0° | 0° | 0° | R | 41-2160-0077 | — |
| | | | | L | 41-2260-0077 | — |



Edgewise Non-Convertible - HG .045\"/>

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|--------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018\"x.025" | 0° | 0° | 0° | R | 41-1160-0084 | — |
| | | | | L | 41-1260-0084 | — |
| .022\"x.028" | 0° | 0° | 0° | R | 41-2160-0084 | — |
| | | | | L | 41-2260-0084 | — |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Maxillary Double Tubes 1st Molar



Edgewise Convertible - HG .045"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | 0° | 6°D | 6°D | R | 41-1160-0081 | — |
| | | | | L | 41-1260-0081 | — |
| .022"x.028" | 0° | 6°D | 6°D | R | 41-2160-0081 | — |
| | | | | L | 41-2260-0081 | — |



Edgewise Non-Convertible - HG .045"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .022"x.028" | 0° | 6°D | 0° | R | 41-2160-0085 | — |
| | | | | L | 41-2260-0085 | — |



Type Planché Convertible - HG .045"/Gingival

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .022"x.028" | -14° | 14°D | 0° | R | 41-2160-0070 | — |
| | | | | L | 41-2260-0070 | — |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Maxillary Triple Tubes 1st Molar



MAXILLARY Triple Tubes 1st Molar
GET HG .045"/Occlusal



| | | | | | | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|--|
| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | Weldable | Bondable | |
| .022"x.028" | -20° | 10°D | 0° | R | 41-2160-0500 | — | |
| | | | | L | 41-2260-0500 | — | |
| .022"x.028" | -20° | 0° | 0° | R | 41-2160-0501 | — | |
| | | | | L | 41-2260-0501 | — | |



Type Roth Convertible - HG .045"/Occlusal

| | | | | | | P/N | |
|-----------------------------|--------|---------------|-------------|-----|--------------|----------|--|
| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | Weldable | Bondable | |
| .018"x.025"/ .018"x.025" | -14° | 14°D | 14°D | R | 41-1160-0010 | — | |
| | | | | L | 41-1260-0010 | — | |
| .022"x.028"/ .018"x.025" | -14° | 14°D | 14°D | R | 41-2160-0010 | — | |
| | | | | L | 41-2260-0010 | — | |



Type Roth Convertible - HG .045"/Occlusal

| | | | | | | P/N | |
|-----------------------------|--------|---------------|-------------|-----|--------------|----------|--|
| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | Weldable | Bondable | |
| .022"x.028"/ .022"x.028" | -14° | 14°D | 14°D | R | 41-2160-0050 | — | |
| | | | | L | 41-2260-0050 | — | |



Type Roth Convertible - HG .045"/Occlusal

| | | | | | | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|--|
| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | Weldable | Bondable | |
| .018"x.025" | -14° | 0° | 0° | R | 41-1160-0051 | — | |
| | | | | L | 41-1260-0051 | — | |
| .022"x.028" | -14° | 0° | 0° | R | 41-2160-0051 | — | |
| | | | | L | 41-2260-0051 | — | |



Type Roth Convertible - HG .051"/Occlusal

| | | | | | | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|--|
| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | Weldable | Bondable | |
| .022"x.028" | -14° | 14°D | 14°D | R | 41-2160-0052 | — | |
| | | | | L | 41-2260-0052 | — | |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Maxillary Triple Tubes 1st Molar



Type Roth Convertible - HG .051"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .022"x.028" | -14° | 0° | 0° | R | 41-2160-0053 | — |
| | | | | L | 41-2260-0053 | — |



Type MBT* Convertible F.E.O. .045"/Occlusal

| Dimensions | Torque | Rotation Slot | Slot aux. | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -14° | 10° | .018"x.025" | R | 41-1160-0111 | — |
| | | | | L | 41-1260-0111 | — |
| .022"x.028" | -14° | 10° | .018"x.025" | R | 41-2160-0111 | — |
| | | | | L | 41-2260-0111 | — |



Type Bio-Progressive Convertible - HG .045"/Gingival

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | 0° | 15°D | 15°D | R | 41-1160-0042 | — |
| | | | | L | 41-1260-0042 | — |



Type Bio-Progressive Convertible - HG .045"/Gingival

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | 0° | 0° | 0° | R | 41-1160-0041 | — |
| | | | | L | 41-1260-0041 | — |
| .022"x.028" | 0° | 0° | 0° | R | 41-2160-0041 | — |
| | | | | L | 41-2260-0041 | — |



Type Bio-Progressive Non-Convertible - HG .045"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | 0° | 0° | 0° | R | 41-1160-0040 | — |
| | | | | L | 41-1260-0040 | — |
| .022"x.028" | 0° | 0° | 0° | R | 41-2160-0040 | — |
| | | | | L | 41-2260-0040 | — |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs


Edgewise Convertible - HG .045"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | 0° | 6°D | 0° | R | 41-1160-0065 | — |
| | | | | L | 41-1260-0065 | — |
| .022"x.028" | 0° | 6°D | 0° | R | 41-2160-0065 | — |
| | | | | L | 41-2260-0065 | — |


Type Polar Convertible - HG .045"/Gingival

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -20° | 10°D | 0° | R | 41-1160-0027 | — |
| | | | | L | 41-1260-0027 | — |


Type Hilgers Convertible - HG .045"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -10° | 15°D | 0° | R | 41-1160-0046 | — |
| | | | | L | 41-1260-0046 | — |
| .022"x.028" | -10° | 15°D | 0° | R | 41-2160-0046 | — |
| | | | | L | 41-2260-0046 | — |


Type Melsen Convertible - HG .045"/Occlusal

| Dimensions | Torque | Rotation Slot | Rotation HG | R/L | P/N | |
|-------------|--------|---------------|-------------|-----|--------------|----------|
| | | | | | Weldable | Bondable |
| .018"x.025" | -10° | 8°D | 0° | R | 41-1160-0068 | — |
| | | | | L | 41-1260-0068 | — |
| .022"x.028" | -10° | 8°D | 0° | R | 41-2160-0068 | — |
| | | | | L | 41-2260-0068 | — |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Mandibular Single Tubes 1st Molar



Convertible



Type Roth

| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -25° | 0° | 4°D | 4.0 | R | 41-1460-0012 | 40-1460-0012 |
| | | | | | L | 41-1360-0012 | 40-1360-0012 |
| .018"x.025" | -30° | 0° | 4°D | 4.0 | R | 41-1460-0062 | 40-1460-0062 |
| | | | | | L | 41-1360-0062 | 40-1360-0062 |
| .022"x.028" | -25° | 0° | 4°D | 4.0 | R | 41-2460-0012 | 40-2460-0012 |
| | | | | | L | 41-2360-0012 | 40-2360-0012 |
| .022"x.028" | -30° | 0° | 4°D | 4.0 | R | 41-2460-0062 | 40-2460-0062 |
| | | | | | L | 41-2360-0062 | 40-2360-0062 |



Type MBT*

| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -20° | 0° | 0° | 3.8 | R | 41-1460-0022 | 40-1460-0119 |
| | | | | | L | 41-1360-0022 | 40-1360-0119 |
| .018"x.025" | -10° | 0° | 0° | 4.0 | R | 41-1460-0076 | 40-1460-0076 |
| | | | | | L | 41-1360-0076 | 40-1360-0076 |
| .022"x.028" | -20° | 0° | 0° | 3.8 | R | 41-2460-0119 | 40-2460-0119 |
| | | | | | L | 41-2360-0119 | 40-2360-0119 |
| .022"x.028" | -10° | 0° | 0° | 4.0 | R | 41-2460-0076 | 40-2460-0076 |
| | | | | | L | 41-2360-0076 | 40-2360-0076 |



Type Alexander

| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|----------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -22° | -6° | 5°D | 3.8 | R | 41-1460-0074 | — |
| | | | | | L | 41-1360-0074 | — |



Type Andrews

| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -25° | 0° | 0° | 3.6 | R | 41-1460-0034 | 40-1460-0034 |
| | | | | | L | 41-1360-0034 | 40-1360-0034 |
| .022"x.028" | -25° | 0° | 0° | 3.6 | R | 41-2460-0034 | 40-2460-0034 |
| | | | | | L | 41-2360-0034 | 40-2360-0034 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

Convertible



Type Planché

| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|----------|
| | | | | | | Weldable | Bondable |
| .022"x.028" | -25° | -6° | 6°D | 4.0 | R | 41-2460-0121 | — |
| | | | | | L | 41-2360-0121 | — |



Edgewise

| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | P/N | |
|-------------|--------|------------|----------|-----------|-----|----------|--------------|
| | | | | | | Weldable | Bondable |
| .022"x.028" | 0° | 0° | 6°D | 4.0 | R | — | 40-2460-0092 |
| | | | | | L | — | 40-2360-0092 |



Edgewise

| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | 0° | 0° | 0° | 4.0 | R | 41-1460-0093 | 40-1460-0093 |
| | | | | | L | 41-1360-0093 | 40-1360-0093 |
| .022"x.028" | 0° | 0° | 0° | 4.0 | R | 41-2460-0093 | 40-2460-0093 |
| | | | | | L | 41-2360-0093 | 40-2360-0093 |



Type Polar

| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | P/N | |
|-------------|--------|------------|----------|-----------|-----|----------|--------------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -22° | 0° | 4°D | 4.0 | R | — | 40-2460-0029 |
| | | | | | L | — | 40-2360-0029 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Mandibular Single Tubes 2nd Molar



Non-Convertible



Type Roth

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -25° | 0° | 6°D | 4.0 | R | 41-1470-0012 | 40-1470-0012 | |
| | | | | | L | 41-1370-0012 | 40-1370-0012 | |
| .022"x.028" | -25° | 0° | 6°D | 4.0 | R | 41-2470-0012 | 40-2470-0012 | |
| | | | | | L | 41-2370-0012 | 40-2370-0012 | |



Type Roth

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -30° | -1° | 4°D | 4.0 | R | 41-1470-0064 | 40-1470-0064 | |
| | | | | | L | 41-1370-0064 | 40-1370-0064 | |
| .022"x.028" | -30° | -1° | 4°D | 4.0 | R | 41-2470-0064 | 40-2470-0064 | |
| | | | | | L | 41-2370-0064 | 40-2370-0064 | |



Type MBT*

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -10° | 0° | 0° | 3.2 | R | 41-1470-0022 | 40-1470-0022 | |
| | | | | | L | 41-1370-0022 | 40-1370-0022 | |
| .022"x.028" | -10° | 0° | 0° | 3.2 | R | 41-2470-0022 | 40-2470-0022 | |
| | | | | | L | 41-2370-0022 | 40-2370-0022 | |



Plan. Type

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|----------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .022"x.028" | -27° | -10° | 6°D | | R | 41-2470-0012 | — | |
| | | | | | L | 41-2370-0012 | — | |



Alexander Type

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|----------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .022"x.028" | -27° | -10° | 6°D | | R | — | 40-2470-0074 | |
| | | | | | L | — | 40-2370-0074 | |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

Non-Convertible



Type Andrews

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|----------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -25° | 0° | 0° | 4.0 | R | — | 40-1470-0034 | |
| | | | | | L | — | 40-1370-0034 | |
| .022"x.028" | -25° | 0° | 0° | 4.0 | R | — | 40-2470-0034 | |
| | | | | | L | — | 40-2370-0034 | |
| .018"x.025" | -25° | 0° | 4°D | 3.0 | R | — | 40-1470-0049 | |
| | | | | | L | — | 40-1370-0049 | |
| .022"x.028" | 25° | 0° | 4°D | 3.0 | R | — | 40-2470-0049 | |
| | | | | | L | — | 40-2370-0049 | |



Type Andrews Wide Entry (Peucat)

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -25° | 0° | 4°D | 3.5 | R | 41-1470-0018 | 40-1470-0018 | |
| | | | | | L | 41-1370-0018 | 40-1370-0018 | |
| .022"x.028" | -25° | 0° | 4°D | 3.5 | R | 41-2470-0018 | 40-2470-0018 | |
| | | | | | L | 41-2370-0018 | 40-2370-0018 | |




Type Polar

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|----------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -25° | 5° | 7°D | 4 | R | — | 40-1470-0029 | |
| | | | | | L | — | 40-1370-0029 | |



Edgewise

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | 0° | 0° | 0° | 4.6 | R | — | 40-1470-0094 | |
| | | | | | L | — | 40-1370-0094 | |
| .022"x.028" | 0° | 0° | 0° | 4.6 | R | 41-2470-0094 | 40-2470-0094 | |
| | | | | | L | 41-2370-0094 | 40-2370-0094 | |



Edgewise

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | 0° | 0° | 0° | 4.0 | R | 41-1470-0095 | 40-1470-0095 | |
| | | | | | L | 41-1370-0095 | 40-1370-0095 | |
| .022"x.028" | 0° | 0° | 0° | 4.0 | R | 41-2470-0095 | 40-2470-0095 | |
| | | | | | L | 41-2370-0095 | 40-2370-0095 | |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Mandibular Single Tubes 2nd Molar



Non-Convertible



Edgewise

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | 0° | 0° | 6°D | 4.0 | R | — | 40-1470-0096 | |
| | | | | | L | — | 40-1370-0096 | |
| .022"x.028" | 0° | 0° | 6°D | 4.0 | R | 41-2470-0096 | 40-2470-0096 | |
| | | | | | L | 41-2370-0096 | 40-2370-0096 | |



Edgewise

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | 0° | 0° | 6°D | 3.2 | R | — | 40-1470-0097 | |
| | | | | | L | — | 40-1370-0097 | |
| .022"x.028" | 0° | 0° | 6°D | 3.2 | R | 41-2470-0096 | 40-2470-0097 | |
| | | | | | L | 41-2370-0096 | 40-2370-0097 | |




Edgewise

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | 0° | 0° | 0° | 3.2 | R | — | 40-1470-0098 | |
| | | | | | L | — | 40-1370-0098 | |
| .022"x.028" | 0° | 0° | 0° | 3.2 | R | 41-2470-0098 | 40-2470-0098 | |
| | | | | | L | 41-2370-0098 | 40-2370-0098 | |



Small

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|----------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Weldable | Bondable | |
| .018"x.025" | -25° | 0° | 4°D | 3.0 | R | — | 40-1470-0048 | |
| | | | | | L | — | 40-1370-0048 | |
| .022"x.028" | -25° | 0° | 4°D | 3.0 | R | — | 40-2470-0048 | |
| | | | | | L | — | 40-2370-0048 | |



Mini Eruption Tube for 2nd Molar

| | | | | | | | P/N | |
|-------------|--------|------------|----------|-----------|-----|--------------|--------------|--|
| Dimensions | Torque | Angulation | Rotation | Length mm | R/L | Without hook | With hook | |
| .018"x.025" | -25° | 0° | 0° | 2.0 | R | 40-1470-0112 | 40-1470-0019 | |
| | | | | | L | 40-1370-0112 | 40-1370-0019 | |
| .022"x.028" | -25° | 0° | 0° | 2.0 | R | 40-2470-0112 | 40-2470-0019 | |
| | | | | | L | 40-2370-0112 | 40-2370-0019 | |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Mandibular Double Tubes 1st Molar



Convertible with Tube Bumper



Type BiDimensional .051"/Gingival

| Dimensions | Torque | Angulation Tube | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|--------|--------------------|------------------|------------------|-----|--------------|----------|
| | | | | | | Weldable | Bondable |
| .022"x.028" | 0° | 0° | 0° | 0° | R | 41-2460-0067 | — |
| | | | | | L | 41-2360-0067 | — |



Type Cetlin .045"/Gingival

| Dimensions | Torque | Angulation Tube | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|--------|--------------------|------------------|------------------|-----|--------------|----------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -25° | 0° | 4°D | 5M | R | 41-1460-0037 | — |
| | | | | | L | 41-1360-0037 | — |
| .022"x.028" | -25° | 0° | 4°D | 5M | R | 41-2460-0037 | — |
| | | | | | L | 41-2360-0037 | — |

Convertible



Type Roth Convertible

| Dimensions | Torque | Angulation Tube | Rotation Slot | Rotation Aux. | R/L | P/N | |
|---------------|--------|--------------------|------------------|------------------|-----|--------------|--------------|
| | | | | | | Weldable | Bondable |
| .018"x.025" - | 25° | 0° | 4°D | 0° | R | 41-1460-0011 | 40-1460-0011 |
| | | | | | L | 41-1360-0011 | 40-1360-0011 |
| .022"x.028" | -25° | 0° | 4°D | 0° | R | 41-2460-0011 | 40-2460-0011 |
| | | | | | L | 41-2360-0011 | 40-2360-001 |



Type Roth Occlusal Hook

| Dimensions | Torque | Angulation Tube | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|--------|--------------------|------------------|------------------|-----|--------------|--------------|
| | | | | | | Weldable | Bondable |
| .022"x.028" | -30° | -1° | 4°D | 0° | R | 41-2460-0058 | 40-2460-0058 |
| | | | | | L | 41-2360-0058 | 40-2360-0058 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs



Mandibular Double Tubes 1st Molar



Convertible



Type Roth Gingival Hook

| Dimensions | Torque | Angulation Tube | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|--------|-----------------|---------------|---------------|-----|--------------|--------------|
| | | | | | | Weldable | Bondable |
| .022"x.028" | -30° | -1° | 4°D | 0° | R | 41-2460-0060 | 40-2460-0060 |
| | | | | | L | 41-2360-0060 | 40-2360-0060 |



Extra Torque

| Dimensions | Torque | Angulation Tube | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|--------|-----------------|---------------|---------------|-----|--------------|--------------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -10° | 0° | 4°D | 0° | R | 41-1460-0070 | — |
| | | | | | L | 41-1360-0070 | — |
| .022"x.028" | -10° | 0° | 4°D | 0° | R | 41-2460-0070 | 40-2460-0070 |
| | | | | | L | 41-2360-0070 | 40-2360-0070 |



Extra Torque Wide Base

| Dimensions | Torque Slot & Aux. | Angulation Tube | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|--------------------|-----------------|---------------|---------------|-----|----------|--------------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -10° | 0° | 4°D | 0° | R | — | 40-1460-0031 |
| | | | | | L | — | 40-1360-0031 |
| .022"x.028" | -10° | 0° | 4°D | 0° | R | — | 40-2460-0031 |
| | | | | | L | — | 40-2360-0031 |



Type Bio-Progressive

| Dimensions | Torque Slot | Torque Aux. | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|-------------|-------------|---------------|---------------|-----|--------------|--------------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | 0° | 0° | 0° | 0° | UNI | 41-1660-0045 | 40-1660-0045 |
| .022"x.028" | 0° | 0° | 0° | 0° | UNI | 41-2660-0045 | 40-2660-0045 |



Type Polar

| Dimensions | Torque Slot | Torque Aux. | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|-------------|-------------|---------------|---------------|-----|--------------|----------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -25° | 0° | 4°D | 0° | R | 41-1460-0028 | — |
| | | | | | L | 41-1360-0028 | — |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

Convertible



Type Hilgers

| Dimensions | Torque Slot | Torque Aux. | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|-------------|-------------|---------------|---------------|-----|--------------|----------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -27° | 0° | 12°D | 0° | R | 41-1460-0047 | — |
| | | | | | L | 41-1360-0047 | — |
| .022"x.028" | -27° | 0° | 12°D | 0° | R | 41-2460-0047 | — |
| | | | | | L | 41-2360-0047 | — |



Type Melsen

| Dimensions | Torque Slot | Torque Aux. | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|-------------|-------------|---------------|---------------|-----|--------------|----------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | -27° | 0° | 8°D | 0° | R | 41-1460-0069 | — |
| | | | | | L | 41-1360-0069 | — |
| .022"x.028" | -27° | 0° | 8°D | 0° | R | 41-2460-0069 | — |
| | | | | | L | 41-2360-0069 | — |

Non- Convertible



Type Roth

| Dimensions | Torque | Angulation Tube | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|--------|-----------------|---------------|---------------|-----|--------------|----------|
| | | | | | | Weldable | Bondable |
| .022"x.028" | -30° | -1° | 4°D | 0° | R | 41-2460-0059 | — |
| | | | | | L | 41-2360-0059 | — |

Type Bio-Progressive



| Dimensions | Torque Slot | Torque Aux. | Rotation Slot | Rotation Aux. | R/L | P/N | |
|-------------|-------------|-------------|---------------|---------------|-----|--------------|--------------|
| | | | | | | Weldable | Bondable |
| .018"x.025" | 0° | 0° | 0° | 0° | UNI | 41-1660-0044 | 40-1660-0044 |
| .022"x.028" | 0° | 0° | 0° | 0° | UNI | 41-2660-0044 | 40-2660-0044 |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs




VIPER Tubes

Buccal Single Non-Convertible

VIPER™ LOW PROFILE BUCCAL TUBES

Innovative technology and sleek, stylish comfort

VIPER™ Low Profile Buccal Tubes Sleek and stylish, G&H®'s Viper™ buccal tube incorporates the latest design technology with a growing industry demand for comfortable, easy-to-use appliances that simplify treatment and reduce valuable chair time.



MAXILLARY 1st & 2nd Molar

| Torque | Offset | Rx | M/D | R/L | ID | Bondable | | | |
|--------|--------|------|--------|-----|----|-----------|-----------|-----------|-----------|
| | | | | | | .018 | .022 | .018 | .022 |
| 0° | 0° | STD | 4.2 mm | R | ◆ | NS3B-821E | NS3B-221E | NS3B-821I | NS3B-221I |
| | | | | L | ◆ | NS3B-822E | NS3B-222E | NS3B-822I | NS3B-222I |
| -10° | 8° | Roth | 4.2 mm | R | ● | NS3B-825E | NS3B-225E | NS3B-825I | NS3B-225I |
| | | | | L | ● | NS3B-826E | NS3B-226E | NS3B-826I | NS3B-226I |
| -14° | 10° | MBT | 4.2 mm | R | ● | NS3B-833E | NS3B-233E | NS3B-833I | NS3B-233I |
| | | | | L | ● | NS3B-834E | NS3B-234E | NS3B-834I | NS3B-234I |

MANDIBULAR 1st & 2nd Molar

| Torque | Offset | Rx | M/D | R/L | ID | Bondable | | | |
|--------|--------|------|--------|-----|----|-----------|-----------|-----------|-----------|
| | | | | | | .018 | .022 | .018 | .022 |
| 0° | 0° | STD | 4.2 mm | R | ◆ | NS3B-822E | NS3B-222E | NS3B-822I | NS3B-222I |
| | | | | L | ◆ | NS3B-821E | NS3B-221E | NS3B-821I | NS3B-221I |
| -25° | 5° | Roth | 4.2 mm | R | ■ | NS3B-827E | NS3B-227E | NS3B-827I | NS3B-227I |
| | | | | L | ■ | NS3B-828E | NS3B-228E | NS3B-828I | NS3B-228I |
| -20° | 0° | MBT | 4.2 mm | R | ▲ | NS3B-809E | NS3B-209E | NS3B-809I | NS3B-209I |
| | | | | L | ▲ | NS3B-810E | NS3B-210E | NS3B-810I | NS3B-210I |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

Evolve Flow Bondable Tubes

Positioning Guide


Occlusal/gingival positioning guides ensure a secure grip and ease of placement

Distal Notch

Evolve Flow features a distal notch for ligation and tieback of auxiliaries

Contoured Base

Anatomically contoured base with occlusal indent ensures accurate fit, whilst the 80 gauge mesh pad ensures superior bond strength



Low profile

The contoured surfaces are low profile and tapered for patient comfort

Low profile hook

Low profile hooks guarantee patient comfort and allow for easy placements of attachments

The Funnelled Entrance

The Evolve Flow funnelled entrance allows for easy wire insertion

| Roth .022 | 1 st Molars | 2 nd Molars |
|-------------|------------------------|------------------------|
| Upper Left | DB07-1350UL | DB07-1371UL |
| Upper Right | DB07-1350UR | DB07-1371UR |
| Lower Left | DB07-1350LL | DB07-1371LL |
| Lower Right | DB07-1350LR | DB07-1371LR |

| MBT .022 | 1 st Molars | 2 nd Molars |
|-------------|------------------------|------------------------|
| Upper Left | DB07-1360UL | DB07-1372UL |
| Upper Right | DB07-1360UR | DB07-1372UR |
| Lower Left | DB07-1360LL | DB07-1372LL |
| Lower Right | DB07-1360LR | DB07-1372LR |

| Roth .018 | 1 st Molars | 2 nd Molars |
|-------------|------------------------|------------------------|
| Upper Left | DB07-1370UL | DB07-1373UL |
| Upper Right | DB07-1370UR | DB07-1373UR |
| Lower Left | DB07-1370LL | DB07-1373LL |
| Lower Right | DB07-1370LR | DB07-1373LR |

| MBT .018 | 1 st Molars | 2 nd Molars |
|-------------|------------------------|------------------------|
| Upper Left | DB07-1362UL | DB07-1374UL |
| Upper Right | DB07-1362UR | DB07-1374UR |
| Lower Left | DB07-1362LL | DB07-1374LL |
| Lower Right | DB07-1362LR | DB07-1374LR |

De tubes zijn verpakt per 10 stuks/ Les tubes sont emballés par 10 pcs

