

Rhodium BRACKETS





EXPERIENCE mini Rhodium



Experience mini RHODIUM Type Roth

MAXILLARY

Tooth	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
U1	12°T	5°A	0°	24-1110-0010	24-1210-0010	24-2110-0010	24-2210-0010
U2	8°T	9°A	0°	24-1120-0010	24-1220-0010	24-2120-0010	24-2220-0010
U3/Hk	-2°T	13°A	4°m	24-113D-0010	24-123D-0010	24-213D-0010	24-223D-0010
U4-5	-7°T	0°A	2°d	24-1140-0010	24-1240-0010	24-2140-0010	24-2240-0010
U4-5/Hk distal	-7°T	0°A	2°d	24-114D-0010	24-124D-0010	24-214D-0010	24-224D-0010

MANDIBULAR

Tooth	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
L1-2	-1°T	2°A	0°	24-1410-0010	24-1310-0010	24-2410-0010	24-2310-0010
L3/Hk	-11°T	7°A	2°m	24-143D-0010	24-133D-0010	24-243D-0010	24-233D-0010
L4	-17°T	-1°A	4°d	24-1440-0010	24-1340-0010	24-2440-0010	24-2340-0010
L4/Hk	-17°T	-1°A	4°d	24-144D-0010	24-134D-0010	24-244D-0010	24-234D-0010
L5	-22°T	-1°A	4°d	24-1450-0010	24-1350-0010	24-2450-0010	24-2350-0010
L5/Hk	-22°T	-1°A	4°d	24-145D-0010	24-135D-0010	24-245D-0010	24-235D-0010

Experience mini RHODIUM Type MBT* *MBT is a trade of 3M Unitek. All rights reserved.

MAXILLARY

Tooth	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
U1	17°	4°	0°	24-1110-0020	24-1210-0020	24-2110-0020	24-2210-0020
U2	10°	8°	0°	24-1120-0020	24-1220-0020	24-2120-0020	24-2220-0020
U3/Hk	-7°	8°	0°	24-113D-0020	24-123D-0020	24-213D-0020	24-223D-0020
U3/Hk/Opt	0°	8°	0°	24-113D-0021	24-123D-0021	24-213D-0021	24-223D-0021
U4-5/Hk	-7°	0°	0°	24-114D-0020	24-124D-0020	24-214D-0020	24-224D-0020
U4-5	-7°	0°	0°	24-1540-0020	24-1540-0020	24-2540-0020	24-2540-0020
U5/Hk/(Opt)	-7°	0°	0°	24-115D-0020	24-125D-0020	24-215D-0020	24-225D-0020
U5/(Opt)	-7°	0°	0°	24-1550-0020	24-1550-0020	24-2550-0020	24-2550-0020

MANDIBULAR

Tooth	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
L1-2	-6°	0°	0°	24-1610-0020	24-1610-0020	24-2610-0020	24-2610-0020
U3/Hk	-6°	3°	0°	24-143D-0020	24-133D-0020	24-243D-0020	24-233D-0020
U3/Hk/Opt	0°	3°	0°	24-143D-0021	24-133D-0021	24-243D-0021	24-233D-0021
L4	-12°	2°	0°	24-1440-0020	24-1340-0020	24-2440-0020	24-2340-0020
L4/Hk	-12°	2°	0°	24-144D-0020	24-134D-0020	24-244D-0020	24-234D-0020
L5	-17°	2°	0°	24-1450-0020	24-1350-0020	24-2450-0020	24-2350-0020
L5/Hk	-17°	2°	0°	24-145D-0020	24-135D-0020	24-245D-0020	24-235D-0020

OPTIONAL BRACKETS

Tooth	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
U4-5	0°	0°	0°	24-1540-0025	24-1540-0025	24-2540-0025	24-2540-0025



EXPERIENCE mini Rhodium



Experience mini RHODIUM, Type Roth .020" x .028"

MAXILLARY

Tooth	Torque	Angulation	Rotation	R	L
U 1 .020" x .028"	12°	5°	0°	24-5110-0010	24-5210-0010
U 2 .020" x .028"	8°	9°	0°	24-5120-0010	24-5220-0010

MANDIBULAR

Tooth	Torque	Angulation	Rotation	R	L
L 1/2 .020" x .028"	-1°	2°	0°	24-5410-0010	24-5310-0010

Experience mini RHODIUM

MAXILLARY



Tooth	Torque	Angulation	Rotation	R	L
U1 .020" x .028"	12°	5°	0°	24-5110-0500	24-5210-0500
U2 .020" x .028"	8°	9°	0°	24-5120-0500	24-5220-0500
U3 .022" x .028"	0°	11°	0°	24-213D-0500	24-223D-0500
U4-5 .022" x .028"	-10°	4°	2°d	24-2140-0500	24-2240-0500
U4-5/Hk .022" x .028"	-10°	4°	2°d	24-214M-0500	24-224M-0500

MANDIBULAR


Tooth	Torque	Angulation	Rotation	R	L
L1-2 .020" x .028"	-1°	2°	0°	24-5410-0500	24-5310-0500
L3/Hk .022" x .028"	-3°	3°	0°	24-243D-0500	24-233D-0500
L4 .022" x .028"	-17°	1°	4°d	24-2440-0500	24-2340-0500
L4/Hk .022" x .028"	-17°	1°	4°d	24-244D-0500	24-234D-0500
L5 .022" x .028"	-22°	1°	4°d	24-2450-0500	24-2350-0500
L5/Hk .022" x .028"	-22°	1°	4°d	24-245D-0500	24-235D-0500

Rhodium brackets zijn verpakt per 10 stuks / Les brackets Rhodium sont emballés par 10 pcs

Experience mini RHODIUM Kits

Prescription	Upper arch	Lower arch	Hook on	.018"	.022"
Experience Mini Rhodium ROTH	5-5	5-5	3	24-1CC3-0010	24-2CC3-0010
Experience Mini Rhodium ROTH	5-5	5-5	3,4,5	24-1CC5-0010	24-2CC5-0010
Experience Mini Rhodium MBT	5-5	5-5	3	24-1CC3-0020	24-2CC3-0020
Experience Mini Rhodium MBT	5-5	5-5	3,4,5	24-1CC5-0020	24-2CC5-0020
G.E.T.	5-5	5-5	3	-	24-5CC3-0500
G.E.T.	5-5	5-5	3,4,5	-	24-5CC5-0500
G.E.T.	5-5	5-5	3	-	24-5AA3-0500
G.E.T.	5-5	5-5	3,4,5	-	24-5AA5-0500





EXPERIENCE

metal Rhodium



Experience RHODIUM, Type Roth

MAXILLARY

Tooth	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
U1	12°	5°	0°	19-1110-0010	19-1210-0010	19-2110-0010	19-2210-0010
U2	8°	9°	0°	19-1120-0010	19-1220-0010	19-2120-0010	19-2220-0010
U3/Hk	-2°	13°	4°m	19-113D-0010	19-123D-0010	19-213D-0010	19-223D-0010
U4-5	-7°	0°	2°d	19-1140-0010	19-1240-0010	19-2140-0010	19-2240-0010
U4-5/Hk d	-7°	0°	2°d	19-114D-0010	19-124D-0010	19-214D-0010	19-224D-0010

MANDIBULAR

Tooth	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
L1-2	-1°	2°	0°	19-1410-0010	19-1310-0010	19-2410-0010	19-2310-0010
L3/Hk	-11°	7°	2°m	19-143D-0010	19-133D-0010	19-243D-0010	19-233D-0010
L4	-17°	-1°	4°d	19-1440-0010	19-1340-0010	19-2440-0010	19-2340-0010
L4/Hk	-17°	-1°	4°d	19-144D-0010	19-134D-0010	19-244D-0010	19-234D-0010
L5	-22°	-1°	4°d	19-1450-0010	19-1350-0010	19-2450-0010	19-2350-0010
L5/Hk	-22°	-1°	4°d	19-145D-0010	19-135D-0010	19-245D-0010	19-235D-0010

Experience RHODIUM, Type MBT*

*MBT is a trade of 3M Unitek. All rights reserved.


MAXILLARY

Tooth	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
U1	17°	4°	0°	19-1110-0020	19-1210-0020	19-2110-0020	19-2210-0020
U2	10°	8°	0°	19-1120-0020	19-1220-0020	19-2120-0020	19-2220-0020
U3/Hk	-7°	8°	0°	19-113D-0020	19-123D-0020	19-213D-0020	19-223D-0020
U3/Hk	0°	8°	0°	19-113D-0021	19-123D-0021	19-213D-0021	19-223D-0021
U4-5	-7°	0°	0°	19-1540-0020	19-1540-0020	19-2540-0020	19-2540-0020
U4-5/Hk d	-7°	0°	0°	19-114D-0020	19-124D-0020	19-214D-0020	19-224D-0020

MANDIBULAR

Tooth	Torque	Angulation	Rotation	.018"		.022"	
				R	L	R	L
L1-2	-6°	0°	0°	19-1610-0020	19-1610-0020	19-2610-0020	19-2610-0020
L3/Hk	-6°	3°	0°	19-143D-0020	19-133D-0020	19-243D-0020	19-233D-0020
L3/Hk	0°	3°	0°	19-143D-0021	19-133D-0021	19-243D-0021	19-233D-0021
L4	-12°	2°	0°	19-1440-0020	19-1340-0020	19-2440-0020	19-2340-0020
L4/Hk	-12°	2°	0°	19-144D-0020	19-134D-0020	19-244D-0020	19-234D-0020
L5	-17°	2°	0°	19-1450-0020	19-1350-0020	19-2450-0020	19-2350-0020
L5/Hk	-17°	2°	0°	19-145D-0020	19-135D-0020	19-245D-0020	19-235D-0020

Rhodium brackets zijn verpakt per 10 stuks / Les brackets Rhodium sont emballés par 10 pcs



EXPERIENCE metal Rhodium



Kits Experience RHODIUM

Prescription	Upper arch	Lower arch	Hook on	.018"	.022"
Experience Rhodium ROTH	5-5	5-5	3	19-1CC3-0010	19-2CC3-0010
Experience Rhodium ROTH	5-5	5-5	3,4,5	19-1CC5-0010	19-2CC5-0010
Experience Rhodium MBT	5-5	5-5	3	19-1CC3-0020	19-2CC3-0020
Experience Rhodium MBT	5-5	5-5	3,4,5	19-1CC5-0020	19-2CC5-0020



What is GET?

GOALS

What are our ideal aesthetic, functional and occlusal treatment goals?

EVIDENCE

What evidence is available to support our decisions?

TREATMENT

What is the best realistic treatment option considering the boundaries of everyday practice and according to sound biomechanical principles?

SPEAKERS



Dr. Stefano Troiani

Private practice of Orthodontics, Siegen, Germany graduated at the University of Aarhus (DK) orthodontic program and received his MSc in Orthodontics in 2003. He has maintained a private practice in Denmark since 2003 and in Germany since 2008. He has worked as Clinical Assistant Professor at the Dept. of Orthodontics of the University of Aarhus from 2005 to 2007. He has been the scientific director of the Alpe Adria Face Master held at the University of Graz, Ljubljana and Zagreb in 2012-2013. He has been consultant orthodontist at the department of maxillo-facial surgery of the hospital of Sønderborg (Denmark) from 2008 to 2009. He is consultant orthodontist at the department of maxillo-facial surgery of the university hospital of Odense (Denmark) from 2008 to present.


[f](#) follow him on facebook Kieferorthopäde Stefano Troiani



Dr. Marc Geserik

Graduated in Dentistry at the University of Ulm (D) and his Dentist training at the University of Basel (CH). In addition, he completed a law degree in medical law.

Since 2005 he runs his private practice in Ulm (D) and since 2006 worldwide pursued his lecturer activity, connected with numerous publications. He spent some times abroad and took place at the University of Kandy (Sri Lanka), the University of Tsurumi (Japan) and the University of Southern California in the USA. He also completed that Further education program Bioesthetic Level I-II and Postgraduate Roth / Williams 2 years program and postgraduate SWING Philosophy under dr. Javier Frenk.



LEGEND medium Rhodium



Legend Medium RHODIUM, Roth Type with vertical slot

MAXILLARY

Tooth	Torque	Angulation	.018"		.022"	
			R	L	R	L
U1	12°T	5°A	39-1110-0034	39-1210-0034	39-2110-0034	39-2210-0034
U2	8°T	9°A	39-1120-0034	39-1220-0034	39-2120-0034	39-2220-0034
U3/Hk	0°T	7°A	39-113D-0040	39-123D-0040	39-213D-0040	39-223D-0040
U4-5	-7°T	0°A	39-1540-0034	39-1540-0034	39-2540-0034	39-2540-0034
UR4/5 HK	-7°T	0°A	39-114M-0034	39-124M-0034	39-214M-0034	39-224M-0034

MANDIBULAR

Tooth	Torque	Angulation	.018"		.022"	
			R	L	R	L
L1-2	0°T	0°A	39-1610-0034	39-1610-0034	39-2610-0034	39-2610-0034
L3/Hk	0°T	5°A	39-143D-0040	39-133D-0040	39-243D-0040	39-233D-0040
L4	-11°T	0°A	39-1640-0034	39-1640-0034	39-2640-0034	39-2640-0034
L5	-17°T	0°A	39-1650-0034	39-1650-0034	39-2650-0034	39-2650-0034
L4/Hk	-11°T	0°A	39-144D-0034	39-134D-0034	39-244D-0034	39-234D-0034
L5/Hk	-17°T	0°A	39-145D-0034	39-135D-0034	39-245D-0034	39-235D-0034

Rhodium brackets zijn verpakt per 10 stuks / Les brackets Rhodium sont emballés par 10 pcs

Legend Medium RHODIUM, MBT Type with vertical slot

*MBT is a trade of 3M Unitek. All rights reserved.

MAXILLARY

Tooth	Torque	Angulation	.018"		.022"	
			R	L	R	L
U1	17°T	5°A	39-1110-0029	39-1210-0029	39-2110-0029	39-2210-0029
U2	8°T	9°A	39-1120-0034	39-1220-0034	39-2120-0034	39-2220-0034
U3/Hk	0°T	7°A	39-113D-0040	39-123D-0040	39-213D-0040	39-223D-0040
U4-5	-7°T	0°A	39-1540-0034	39-1540-0034	39-2540-0034	39-2540-0034
UR4/5 HK	-7°T	0°A	39-114M-0034	39-124M-0034	39-214M-0034	39-224M-0034

MANDIBULAR

Tooth	Torque	Angulation	.018"		.022"	
			R	L	R	L
L1-2	-5°T	0°A	39-1610-0011	39-1610-0011	39-2610-0036	39-2610-0036
L3/Hk	0°T	5°A	39-143D-0040	39-133D-0040	39-243D-0040	39-233D-0040
L4	-11°T	0°A	39-1640-0034	39-1640-0034	39-2640-0034	39-2640-0034
L5	-17°T	0°A	39-1650-0034	39-1650-0034	39-2650-0034	39-2650-0034
L4/Hk	-11°T	0°A	39-144D-0034	39-134D-0034	39-244D-0034	39-234D-0034
L5/Hk	-17°T	0°A	39-145D-0034	39-135D-0034	39-245D-0034	39-235D-0034

Rhodium brackets zijn verpakt per 10 stuks / Les brackets Rhodium sont emballés par 10 pcs

LEGEND medium

Rhodium



Legend Medium RHODIUM, Ricketts Type with vertical slot

MAXILLARY

Tooth	Torque	Angulation	.018"		.022"	
			R	L	R	L
U1	22°T	5°A	39-1110-0028	39-1210-0028	39-2110-0028	39-2210-0028
U2	14°T	9°A	39-1120-0028	39-1220-0028	39-2120-0028	39-2220-0028
U3/Hk	7°T	5°A	39-113D-0028	39-123D-0028	39-213D-0028	39-223D-0028
U4-5	0°T	0°A	39-1540-0040	39-1540-0040	39-2540-0040	39-2540-0040

MANDIBULAR

Tooth	Torque	Angulation	.018"		.022"	
			R	L	R	L
L1-2	0°T	0°A	39-1610-0034	39-1610-0034	39-2610-0034	39-2610-0034
L3/Hk	7°T	5°A	39-123D-0028	39-113D-0028	39-223D-0028	39-213D-0028
L4-5	0°T	0°A	39-1640-0040	39-1640-0040	39-2640-0040	39-2640-0040

Rhodium brackets zijn verpakt per 10 stuks / Les brackets Rhodium sont emballés par 10 pcs



LEGEND mini Rhodium



Legend Mini RHODIUM, Roth Type MAXILLARY

Tooth	Torque	Ang.	Rotation	.018"		.022"	
				R	L	R	L
U1	12°	5°	0°	38-1110-0010	38-1210-0010	38-2110-0010	38-2210-0010
U2	8°	9°	0°	38-1120-0010	38-1220-0010	38-2120-0010	38-2220-0010
U3/Hk	2°	10°	4°m	38-113D-0010	38-123D-0010	38-213D-0010	38-223D-0010
U4	-7°	0°	2°d	38-1140-0010	38-1240-0010	38-2140-0010	38-2240-0010
U4/Hk	-7°	0°	2°d	38-114M-0010	38-124M-0010	38-214M-0010	38-224M-0010
U5	-7°	0°	2°d	38-1150-0010	38-1250-0010	38-2150-0010	38-2250-0010
U5/Hk	-7°	0°	2°d	38-115M-0010	38-125M-0010	38-215M-0010	38-225M-0010

MANDIBULAR

Tooth	Torque	Ang.	Rotation	.018"		.022"	
				R	L	R	L
L1-2	-1°	0°	0°	38-1610-0010	38-1610-0010	38-2610-0010	38-2610-0010
L3/Hk	-11°	7°	2°m	38-143D-0010	38-133D-0010	38-243D-0010	38-233D-0010
L4	-17°	0°	4°d	38-1440-0010	38-1340-0010	38-2440-0010	38-2340-0010
L4/Hk	-17°	0°	4°d	38-144D-0010	38-134D-0010	38-244D-0010	38-234D-0010
L5	-22°	0°	4°d	38-1450-0010	38-1350-0010	38-2450-0010	38-2350-0010
L5/Hk	-22°	0°	4°d	38-145D-0010	38-135D-0010	38-245D-0010	38-235D-0010

Legend Mini RHODIUM, MBT *MBT is a trade of 3M Unitek. All rights reserved.

MAXILLARY

Tooth	Torque	Ang.	.018"		.022"	
			R	L	R	L
U1	17°	4°	38-1110-0020	38-1210-0020	38-2110-0020	38-2210-0020
U2	10°	8°	38-1120-0020	38-1220-0020	38-2120-0020	38-2220-0020
U3	-7°	8°	38-1130-0020	38-1230-0020	38-2130-0020	38-2230-0020
U3/Hk	-7°	8°	38-113D-0020	38-123D-0020	38-213D-0020	38-223D-0020
U4	-7°	0°	38-1540-0020	38-1540-0020	38-2540-0020	38-2540-0020
U4/Hk	-7°	0°	38-114D-0020	38-124D-0020	38-214D-0020	38-224D-0020
U5	-7°	0°	38-1550-0020	38-1550-0020	38-2550-0020	38-2550-0020
U5/Hk	-7°	0°	38-115D-0020	38-125D-0020	38-215D-0020	38-225D-0020

MANDIBULAR

Tooth	Torque	Ang.	.018"		.022"	
			R	L	R	L
L 1-2	-6°	0°	38-1610-0020	38-1610-0020	38-2610-0020	38-2610-0020
L 3	-6°	3°	38-1430-0020	38-1330-0020	38-2430-0020	38-2330-0020
L 3/Hk	-6°	3°	38-143D-0020	38-133D-0020	38-243D-0020	38-233D-0020
L 4	-12°	2°	38-1440-0020	38-1340-0020	38-2440-0020	38-2340-0020
L 4/Hk	-12°	2°	38-144D-0020	38-134D-0020	38-244D-0020	38-234D-0020
L 5	-17°	2°	38-1450-0020	38-1350-0020	38-2450-0020	38-2350-0020
L 5/Hk	-17°	2°	38-145D-0020	38-135D-0020	38-245D-0020	38-235D-0020

Rhodium brackets zijn verpakt per 10 stuks / Les brackets Rhodium sont emballés par 10 pcs

Kits Legend Mini RHODIUM

Hook on		.018"	.022"
Mini MBT Type	3	38-1CC3-0020	38-2CC3-0020
Mini MBT Type	3,4,5	38-1CC5-0020	38-2CC5-0020
Mini Roth Type	3	38-1CC3-0010	38-2CC3-0010
Mini Roth Type	3,4,5	38-1CC5-0010	38-2CC5-0010

TYPODONTS

Rhodium



TYPODONT Conventional Rhodium brackets

Legend Mini

	Ref.
With hinge 7-7	55-3811-0000

Legend Medium

	Ref.
With hinge 7-7	55-3911-0000

TYPODONT Self-Ligating Rhodium brackets

EXPERIENCE RHODIUM

	Ref.
With hinge 7-7	55-1911-0000



50-2130-H000

The universel Experience instrument's (handle only)

50-2130-B000

Additional blades (Verpakt per 2 / Emballé par 2)



50-0032-0000

For Experience Mini metal/rhodium



50-0030-0000

For Experience metal/ceramic and rhodium



50-2030-0000

For Experience metal/rhodium

