



Keramische/Céramiques BRACKETS



Chic. Type Roth*



CHIC. Type Roth MAXILLARY

Tooth	Torque	Angulation	P/N			
			.018"		.022"	
			R	L	R	L
U 1	12°	5°	14-1110-0010	14-1210-0010	14-2110-0010	14-2210-0010
U 2	8°	9°	14-1120-0010	14-1220-0010	14-2120-0010	14-2220-0010
U 3	-2°	13°	14-1130-0010	14-1230-0010	14-2130-0010	14-2230-0010
U 3/Hk	-2°	13°	14-113D-0010	14-123D-0010	14-213D-0010	14-223D-0010
U 4-5	-7°	0°	14-1540-0010	14-1540-0010	14-2540-0010	14-2540-0010
U 4-5/Hk	-7°	0°	14-114D-0010	14-124D-0010	14-214D-0010	14-224D-0010

MANDIBULAR

Tooth	Torque	Angulation	P/N			
			.018"		.022"	
			R	L	R	L
L 1-2	0°	0°	14-1610-0010	14-1610-0010	14-2610-0010	14-2610-0010
L 3	-11°	7°	14-1430-0010	14-1330-0010	14-2430-0010	14-2330-0010
L 3/Hk	-11°	7°	14-143D-0010	14-133D-0010	14-243D-0010	14-233D-0010
L 4	-17°	0°	14-1640-0010	14-1640-0010	14-2640-0010	14-2640-0010
L 4/Hk	-17°	0°	14-144D-0010	14-134D-0010	14-244D-0010	14-234D-0010
L 5	-22°	0°	14-1650-0010	14-1650-0010	14-2650-0010	14-2650-0010
L 5/Hk	-22°	0°	14-145D-0010	14-135D-0010	14-245D-0010	14-235D-0010

Keramische brackets zijn verpakt per 5 stuks /
Les brackets céramiques sont emballés par 5 pcs

GC Orthodontics version of the **MBT™ System** prescription
(McLaughlin, Bennet or Trevisi) and is not a copy of existing system.
No endorsement by the Drs McLaughlin, Bennett or Trevisi is implied.

GC Orthodontics version of the **Roth System** prescription
and is not a copy of existing system.
No endorsement by Dr Roth is implied.



Chic. Type MBT*



CHIC. Type MBT*

MAXILLARY

Tooth	I/Omm	Torque	Angulation	P/N		P/N	
				.018"		.022"	
				R	L	R	L
U 1	0.80	17°	4°	14-1110-0020	14-1210-0020	14-2110-0020	14-2210-0020
U 2	1.32	10°	8°	14-1120-0020	14-1220-0020	14-2120-0020	14-2220-0020
U 3/Hk	0.73	-7°	8°	14-113D-0020	14-123D-0020	14-213D-0020	14-223D-0020
U 3/Hk	0.73	0°	8°	14-113D-0021	14-123D-0021	14-213D-0021	14-223D-0021
U 4-5	0.77	-7°	0°	14-1540-0020	14-1540-0020	14-2540-0020	14-2540-0020
U 5/(Optional)	1.07	-7°	0°	14-1550-0020	14-1550-0020	14-2550-0020	14-2550-0020
U 4-5/Hk	0.77	-7°	0°	14-114D-0020	14-124D-0020	14-214D-0020	14-224D-0020

MANDIBULAR

Tooth	I/Omm	Torque	Angulation	P/N		P/N	
				.018"		.022"	
				R	L	R	L
L 1-2	1.26	-6°	0°	14-1610-0020	14-1610-0020	14-2610-0020	14-2610-0020
L 3/Hk	0.83	-6°	3°	14-143D-0020	14-133D-0020	14-243D-0020	14-233D-0020
L 3/Hk	0.83	0°	3°	14-143D-0021	14-133D-0021	14-243D-0021	14-233D-0021
L 4	0.95	-12°	0°	14-1640-0020	14-1640-0020	14-2640-0020	14-2640-0020
L 4/Hk	0.95	-12°	0°	14-144D-0020	14-134D-0020	14-244D-0020	14-234D-0020
L 5	0.95	-17°	0°	14-1650-0020	14-1650-0020	14-2650-0020	14-2650-0020
L 5/Hk	0.95	-17°	0°	14-145D-0020	14-135D-0020	14-245D-0020	14-235D-0020

Keramische brackets zijn verpakt per 5 stuks / Les brackets céramiques sont emballés par 5 pcs

CHIC. KITS

Prescription	Upper	Lower	Hook on	P/N	P/N
				.018"	.022"
Roth	3-3	—	3	14-1E03-0010	14-2E03-0010
Roth	3-3	3-3	3	14-1EE3-0010	14-2EE3-0010
Roth	5-5	—	3	14-1C03-0010	14-2C03-0010
Roth	5-5	3-3	3	14-1CE3-0010	14-2CE3-0010
Roth	5-5	5-5	3	14-1CC3-0010	14-2CC3-0010
Roth	5-5	3-3	3,4,5	14-1CE5-0010	14-2CE5-0010
Roth	5-5	5-5	3,4,5	14-1CC5-0010	14-2CC5-0010
Type MBT	3-3	—	3	—	14-2E03-0020
Type MBT	5-5	—	3	14-1C03-0020	14-2C03-0020
Type MBT	5-5	3-3	3	—	14-2CE3-0020
Type MBT	5-5	5-5	3	14-1CC3-0020	14-2CC3-0020
Type MBT	5-5	3-3	3,4,5	—	14-2CE5-0020
Type MBT	5-5	5-5	3,4,5	14-1CC5-0020	14-2CC5-0020



EXPERIENCE ceramic Type Roth*



EXPERIENCE CERAMIC Type Roth

MAXILLARY

Tooth	I/O mm	Torque	Angulation	Rotation	P/N			
					.018"		.022"	
					R	L	R	L
U 1	.81	12°	5°	0°	11-1110-0010	11-1210-0010	11-2110-0010	11-2210-0010
U 2	1.07	8°	9°	0°	11-1120-0010	11-1220-0010	11-2120-0010	11-2220-0010
U 3	.46	-2°	13°	4°m	11-1130-0010	11-1230-0010	11-2130-0010	11-2230-0010
U 3/Hk	.46	-2°	13°	4°m	11-113D-0010	11-123D-0010	11-213D-0010	11-223D-0010
U 4-5	.46	-7°	0°	2°d	11-1140-0010	11-1240-0010	11-2140-0010	11-2240-0010
U 4-5/Hk	.46	-7°	0°	2°d	11-114D-0010	11-124D-0010	11-214D-0010	11-224D-0010

MANDIBULAR

Tooth	I/O mm	Torque	Angulation	Rotation	P/N			
					.018"		.022"	
					R	L	R	L
L 1-2	1.17	-1°	2°	0°	11-1410-0010	11-1310-0010	11-2410-0010	11-2310-0010
L 3	.56	-11°	7°	2°m	11-1430-0010	11-1330-0010	11-2430-0010	11-2330-0010
L 3/Hk	.56	-11°	7°	2°m	11-143D-0010	11-133D-0010	11-243D-0010	11-233D-0010
L 4(zirconium)	.38	-17°	0°	4°d	11-1440-Z010	11-1340-Z010	11-2440-Z010	11-2340-Z010
L 4/Hk/ (zyrc.)	.38	-17°	0°	4°d	11-144D-Z010	11-134D-Z010	11-244D-Z010	11-234D-Z010
L 5(zirconium)	.28	-22°	0°	4°d	11-1450-Z010	11-1350-Z010	11-2450-Z010	11-2350-Z010
L 5/Hk/(zyrc.)	.28	-22°	0°	4°d	11-145D-Z010	11-135D-Z010	11-245D-Z010	11-235D-Z010

EXPERIENCE CERAMIC Type Roth .020"x .028"

Features an original slot size on upper and lower incisors: .020" x .028". Used in conjunction with .020" x .025" stainless steel archwires, the system controls the anterior torque and brings more rigidity to the buccal segments during space closure.

MAXILLARY

Tooth	I/O mm	Torque	Angulation	Rotation	P/N	
					R	L
U 1 .020" x .028"	.81	12°	5°	0°	11-5110-0010	11-5210-0010
U 2 .020" x .028"	1.07	8°	9°	0°	11-5120-0010	11-5220-0010
U 3 .022" x .028"/Hk	.46	0°	8°	2°m	11-213D-0030	11-223D-0030
U 4/5 .022" x .028"	.46	-7°	0°	2°d	11-2140-0010	11-2240-0010
U 4/5 .022" x .028"/Hk	.46	-7°	0°	2°d	11-214D-0010	11-224D-0010

MANDIBULAR

Tooth	I/O mm	Torque	Angulation	Rotation	P/N	
					R	L
L 1/2 .020" x .028"	1.17	0°	2°	0°	11-5410-0010	11-5310-0010
L 3 .022" x .028"/Hk	.56	-3°	7°	2°m	11-243D-0030	11-233D-0030
L 4* .022" x .028"/Hk/(zirc.)	.38	-17°	0°	4°d	11-244D-Z010	11-234D-Z010
L 5* .022" x .028"/Hk/(zirc.)	.28	-22°	0°	4°d	11-245D-Z010	11-235D-Z010

Keramische brackets zijn verpakt per 5 stuks /
Les brackets céramiques sont emballés par 5 pcs

Keramische/Céramiques BRACKETS



EXPERIENCE ceramic

Extra Torque



EXPERIENCE CERAMIC Extra Torque

This prescription is based on the use of less than full size archwires due to increased torque on the maxillary incisors. Increased buccal root torque in the upper bicuspids prevents any lingual cusp interference.

MAXILLARY

Tooth	I/O mm	Torque	Angulation	Rotation	P/N			
					.018"		.022"	
					R	L	R	L
U 1	.81	18°	5°	0°	11-1110-0030	11-1210-0030	11-2110-0030	11-2210-0030
U 2	1.07	10°	8°	0°	11-1120-0030	11-1220-0030	11-2120-0030	11-2220-0030
U 3	.46	0°	8°	2°m	11-1130-0030	11-1230-0030	11-2130-0030	11-2230-0030
U 3/Hk	.46	0°	8°	2°m	11-113D-0030	11-123D-0030	11-213D-0030	11-223D-0030
U 4-5	.46	-10°	4°	2°d	11-1140-0030	11-1240-0030	11-2140-0030	11-2240-0030
U 4-5/Hk	.46	-10°	4°	2°d	11-114D-0030	11-124D-0030	11-214D-0030	11-224D-0030

MANDIBULAR

Tooth	I/O mm	Torque	Angulation	Rotation	P/N			
					.018"		.022"	
					R	L	R	L
L 1-2	1.17	1°	0°	0°	11-1610-0030	11-1610-0030	11-2610-0030	11-2610-0030
L 3	.56	-3°	7°	2°m	11-1430-0030	11-1330-0030	11-2430-0030	11-2330-0030
L 3/Hk	.56	-3°	7°	2°m	11-143D-0030	11-133D-0030	11-243D-0030	11-233D-0030
L 4*/(Exp M)	.33	-7°	0°	2°d	20-1440-0030	20-1340-0030	20-2440-0030	20-2340-0030
L 4*/Hk/(Exp M)	.33	-7°	0°	2°d	20-144D-0030	20-134D-0030	20-244D-0030	20-234D-0030
L 5*/(Exp M)	.33	-7°	0°	2°d	20-1450-0030	20-1350-0030	20-2450-0030	20-2350-0030
L 5*/Hk/(Exp M)	.33	-7°	0°	2°d	20-145D-0030	20-135D-0030	20-245D-0030	20-235D-0030

* Since ceramic material is much harder than enamel, we recommend using EXPERIENCE metal brackets on mandibular bicuspids (L4, L5)

Keramische brackets zijn verpakt per 5 stuks /
Les brackets céramiques sont emballés par 5 pcs





EXPERIENCE ceramic

Type MBT*



EXPERIENCE CERAMIC Type MBT*

MAXILLARY

Tooth	I/O mm	Torque	Angulation	Rotation	P/N			
					.018"		.022"	
					R	L	R	L
U 1	.81	17°	4°	0°	11-1110-0020	11-1210-0020	11-2110-0020	11-2210-0020
U 2	1.07	10°	8°	0°	11-1120-0020	11-1220-0020	11-2120-0020	11-2220-0020
U 3/Hk	.46	-7°	8°	0°	11-113D-0020	11-123D-0020	11-213D-0020	11-223D-0020
U 3/Hk	.46	0°	8°	0°	11-113D-0021	11-123D-0021	11-213D-0021	11-223D-0021
U 4-5	.46	-7°	0°	0°	11-1540-0020	11-1540-0020	11-2540-0020	11-2540-0020
U 4-5/Hk	.46	-7°	0°	0°	11-114D-0020	11-124D-0020	11-214D-0020	11-224D-0020
U 5/(Optional)	.96	-7°	0°	0°	11-1550-0020	11-1550-0020	11-2550-0020	11-2550-0020
U 5/Hk/(Opt)	.96	-7°	0°	0°	11-115D-0020	11-125D-0020	11-215D-0020	11-225D-0020

MANDIBULAR

Tooth	I/O mm	Torque	Angulation	Rotation	P/N			
					.018"		.022"	
					R	L	R	L
L 1-2	1.17	-6°	0°	0°	11-1610-0020	11-1610-0020	11-2610-0020	11-2610-0020
L 3/Hk	.56	-6°	3°	0°	11-143D-0020	11-133D-0020	11-243D-0020	11-233D-0020
L 3/Hk	.56	0°	3°	0°	11-143D-0021	11-133D-0021	11-243D-0021	11-233D-0021
L 4	.33	-12°	0°	0°	11-1440-0020	11-1340-0020	11-2440-0020	11-2340-0020
L 4/Hk	.33	-12°	0°	0°	11-144D-0020	11-134D-0020	11-244D-0020	11-234D-0020
L 5	.33	-17°	0°	0°	11-1450-0020	11-1350-0020	11-2450-0020	11-2350-0020
L 5/Hk	.33	-17°	0°	0°	11-145D-0020	11-135D-0020	11-245D-0020	11-235D-0020

Keramische brackets zijn verpakt per 5 stuks /
Les brackets céramiques sont emballés par 5 pcs



Keramische/Céramiques BRACKETS



EXPERIENCE ceramic Kits



EXPERIENCE ceramic Kits

				P/N	
Prescription	Upper arch	Lower arch	With Hook	.020" x .028"	
Type Roth .020" x .028"	5-5	3-3	3	11-5CE3-0010	
Type Roth .020" x .028"	5-5	—	3	11-5C03-0010	
				P/N	
Prescription	Upper arch	Lower arch	With Hook	.018"	.022"
Type MBT*	3-3	—	3	11-1E03-0020	11-2E03-0020
Type MBT*	3-3	3-3	3	11-1EE3-0020	11-2EE3-0020
Type MBT*	5-5	—	3	11-1C03-0020	11-2C03-0020
Type MBT*	5-5	—	3,4,5	11-1C05-0020	11-2C05-0020
Type MBT*	5-5	3-3	3	11-1CE3-0020	11-2CE3-0020
Type MBT*	5-5	3-3	3,4,5	11-1CE5-0020	11-2CE5-0020
Type MBT*	5-5	5-5	3	11-1CC3-0020	11-2CC3-0020
Type MBT*	5-5	5-5	3,4,5	11-1CC5-0020	11-2CC5-0020
				P/N	
Prescription	Upper arch	Lower arch	With Hook	.018"	.022"
Type Roth	3-3	—	3	11-1E03-0010	11-2E03-0010
Type Roth	3-3	3-3	3	11-1EE3-0010	11-2EE3-0010
Type Roth	5-5	—	3	11-1C03-0010	11-2C03-0010
Type Roth	5-5	—	3,4,5	11-1C05-0010	11-2C05-0010
Type Roth	5-5	3-3	3	11-1CE3-0010	11-2CE3-0010
Type Roth	5-5	3-3	3,4,5	11-1CE5-0010	11-2CE5-0010
Type Roth (L 4-5 zirconium)	5-5	5-5	3	11-1CC3-M010	11-2CC3-M010
Type Roth (L 4-5 zirconium)	5-5	5-5	3,4,5	11-1CC5-M010	11-2CC5-M010
				P/N	
Prescription	Upper arch	Lower arch	With Hook	.018"	.022"
Extra Torque	3-3	—	3	11-1E03-0030	11-2E03-0030
Extra Torque	3-3	3-3	3	11-1EE3-0030	11-2EE3-0030
Extra Torque	5-5	—	3	11-1C03-0030	11-2C03-0030
Extra Torque	5-5	—	3,4,5	11-1C05-0030	11-2C05-0030
Extra Torque	5-5	3-3	3	11-1CE3-0030	11-2CE3-0030
Extra Torque	5-5	3-3	3,4,5	11-1CE5-0030	11-2CE5-0030
Prescription	Upper arch	Lower arch	With Hook	.022"	
CCO	5-5	-	3,4,5	11-2C05-0007	
CCO	5-5	3-3	3,4,5	11-2CE5-0007	



EXPERIENCE ceramic CCO

EXPERIENCE CERAMIC CCO MAXILLARY

Tooth	I/O mm	Torque	Angulation	P/N .022" x.028"	
				R	L
U 1	.81	12°	5°	11-2110-0007	11-2210-0007
U 2	1.07	10°	9°	11-2120-0007	11-2220-0007
U 3/Hk	.46	-7°	10°	11-213D-0007	11-223D-0007
U 4-5	.46	-9°	0°	11-2540-0007	11-2540-0007
U 4-5/Hk	.46	-9°	0°	11-214D-0007	11-224D-0007

MANDIBULAR

Tooth	I/O mm	Torque	Angulation	P/N .022" x.028"	
				R	L
L 1/2	1.17	-6°	0°	11-2610-0007	11-2610-0007
L 3/Hk	.56	-8°	3°	11-243D-0007	11-233D-0007
L 4	.38	-12°	0°	11-2440-0020	11-2340-0020
L 4/Hk	.38	-12°	0°	11-244D-0020	11-234D-0020
L 5	.38	-17°	0°	11-2450-0020	11-2350-0020
L 5/Hk	.38	-17°	0°	11-245D-0020	11-235D-0020



EXPERIENCE ceramic Optional



OPTIONAL BRACKETS

Tooth	Torque	Angulation	Rotation	P/N			
				.018"	.022"	.022"	.022"
U 3/Hk	7°	5°	0°	11-113D-0027	11-123D-0027	11-213D-0027	11-223D-0027
U 4 Univ	0°	0°	0°	11-1540-0025	11-1540-0025	11-2540-0025	11-2540-0025

Keramische brackets zijn verpakt per 5 stuks / Les brackets céramiques sont emballés par 5 pcs

Instruments for working with EXPERIENCE brackets:



C-Tool: Features both a double forked end to facilitate wire insertion and an end with a hollow ball to safely open the bracket from the gingival side.

P/N. **50-1030-0000**



Open Tool: Is a double ended instrument that features a hollow ball, which aides in opening EXPERIENCE brackets from either the buccal (hook) or gingival (safety) side.

P/N. **50-0030-0000**



20%
smaller than
the Experience
Ceramic

EXPERIENCE mini ceramic Type MBT* en Roth*



EXPERIENCE mini Ceramic Type MBT

MAXILLARY

Tooth	Torque	Angulation	.018		.022	
			R	L	R	L
U1	17°T	4°A	12-1110-0020	12-1210-0020	12-2110-0020	12-2210-0020
U2	10°T	8°A	12-1120-0020	12-1220-0020	12-2120-0020	12-2220-0020
U3 with hook	0°T	8°A	12-113D-0021	12-123D-0021	12-213D-0021	12-223D-0021
U4 univ	-7°T	0°A	12-1540-0020	12-1540-0020	12-2540-0020	12-2540-0020
U4&5 with hook	-7°T	0°A	12-114D-0020	12-124D-0020	12-214D-0020	12-224D-0020

MANDIBULAR

Tooth	Torque	Angulation	.018		.022	
			R	L	R	L
L1-2	-6°T	0°A	12-1610-0020	12-1610-0020	12-2610-0020	12-2610-0020
L3 with hook	0°T	3°A	12-143D-0021	12-133D-0021	12-243D-0021	12-233D-0021
L4	-12°T	2°A	12-1440-0020	12-1340-0020	12-2440-0020	12-2340-0020
L4 with hook	-12°T	2°A	12-144D-0020	12-134D-0020	12-244D-0020	12-234D-0020
L5	-17°T	2°A	12-1450-0020	12-1350-0020	12-2450-0020	12-2350-0020
L5 with hook	-17°T	2°A	12-145D-0020	12-135D-0020	12-245D-0020	12-235D-0020

Keramische brackets zijn verpakt per 5 stuks / Les brackets céramiques sont emballés par 5 pcs

EXPERIENCE mini Ceramic Type Roth

MAXILLARY

Tooth	I/O mm	Torque	Ang.	Rotation	.018		.022	
					R	L	R	L
U1	.81	12°T	5°A	0°	12-1110-0010	12-1210-0010	12-2110-0010	12-2210-0010
U2	1.07	8°T	9°A	0°	12-1120-0010	12-1220-0010	12-2120-0010	12-2220-0010
U3/hook	.46	-2°T	13°A	4°m	12-113D-0010	12-123D-0010	12-213D-0010	12-223D-0010
U4-5/hook	.46	-7°T	0°A	2°d	12-114D-0010	12-124D-0010	12-214D-0010	12-224D-0010

MANDIBULAR

Tooth	I/O mm	Torque	Ang.	Rotation	.018		.022	
					R	L	R	L
L1-2	1.17	-1°T	2°A	0°	12-1410-0010	12-1310-0010	12-2410-0010	12-2310-0010
L3/hook	.56	-11°T	7°A	2°m	12-143D-0010	12-133D-0010	12-243D-0010	12-233D-0010
L4/hook	.38	-17°T	0°A	4°d	12-144D-0010	12-134D-0010	12-244D-0010	12-234D-0010
L5/hook	.28	-22°T	0°A	4°d	12-145D-0010	12-135D-0010	12-245D-0010	12-235D-0010

KITS mini Ceramic

Prescription	Upper	Lower	Hook on	.018	.022
Type Roth	5-5	5-5	3,4,5	12-1CC5-0010	12-2CC5-0010
Type MBT	5-5	5-5	3,4,5	12-1CC5-0020	12-2CC5-0020



AGILITY ceramic Type Roth*

AGILITY CERAMIC Type Roth Sold in packages of 5

MAXILLARY Roth¹ Rx

Tooth	Torque	Angle	Rotation	M/D	R/L	ID	.018	.022
Central	12°	5°	0°	3.9 mm	R	●	127-U1-8R-V	127-U1-2R-V
					L	●	127-U1-8L-V	127-U1-2L-V
Lateral	8°	9°	0°	3.3 mm	R	○	127-U2-8R-V	127-U2-2R-V
					L	○	127-U2-8L-V	127-U2-2L-V
Cuspid (without Hook)	0°	11°	0°	3.4 mm	R	●	127-U3-8R-V	127-U3-2R-V
					L	●	127-U3-8L-V	127-U3-2L-V
Cuspid Hk (with Hook)	0°	11°	0°	3.4 mm	R	●	127-U3-8RK-V	127-U3-2RK-V
					L	●	127-U3-8LK-V	127-U3-2LK-V
1 st & 2 nd Bicuspid	-7°	0°	0°	3.4 mm	U	●●	127-U4-8U-V	127-U4-2U-V
1 st & 2 nd Bicuspid Hk	-7°	0°	0°	3.4 mm	R	●	127-U4-8RK-V	127-U4-2RK-V
					L	●	127-U4-8LK-V	127-U4-2LK-V

MANDIBULAR Roth¹ Rx

Tooth	Torque	Angle	Rotation	M/D	R/L	ID	.018	.022
Anterior	0°	0°	0°	2.6 mm	U	●●	127-L1-8U-V	127-L1-2U-V
Cuspid (without Hook)	-11°	5°	0°	3.4 mm	R	●	127-L3-8R-V	127-L3-2R-V
					L	●	127-L3-8L-V	127-L3-2L-V
Cuspid Hk (with Hook)	-11°	5°	0°	3.4 mm	R	●	127-L3-8RK-V	127-L3-2RK-V
					L	●	127-L3-8LK-V	127-L3-2LK-V
1 st Bicuspid	-17°	0°	0°	3.4 mm	U	●●	127-L4-8U-V	127-L4-2U-V
1 st Bicuspid Hk	-17°	0°	0°	3.4 mm	R	●	127-L4-8RK-V	127-L4-2RK-V
					L	●	127-L4-8LK-V	127-L4-2LK-V
2 nd Bicuspid	-22°	0°	0°	3.4 mm	U	●●	127-L5-8U-V	127-L5-2U-V
2 nd Bicuspid Hk	-22°	0°	0°	3.4 mm	R	●	127-L5-8RK-V	127-L5-2RK-V
					L	●	127-L5-8LK-V	127-L5-2LK-V

BRACKET KITS Roth¹ Rx


20 brackets per kit	.018	.022
Upper/Lower 5x5 Hook on 3	19K-3-18-01	19K-3-22-01
Upper/Lower 5x5 Hook on 3, 4, & 5	19K-345-18-01	19K-345-22-01

Brackets are packaged with 1 case oriented in a kit box.

Brackets are also available in kits of 5 or 10 cases oriented in a carousel dispenser.



Agility™ Twin Bracket
kit, 20 brackets
oriented in a kit box



**MBT is a trademark of 3M Unitek Corporation.

G&H® prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement by the doctor.



AGILITY ceramic Type MBT*

AGILITY CERAMIC Type MBT Sold in packages of 5

MAXILLARY MBT² Rx

Tooth	Torque	Angle	Rotation	M/D	R/L	ID	.018	.022
Central	17°	4°	0°	3.9 mm	R	●	128-U1-8R-V	128-U1-2R-V
					L	●	128-U1-8L-V 1	28-U1-2L-V
Lateral	10°	8°	0°	3.3 mm	R	○	128-U2-8R-V	128-U2-2R-V
					L	○	128-U2-8L-V	128-U2-2L-V
Cuspid (without Hook)	-7°	8°	0°	3.4 mm	R	●	128-U3-8R-V	128-U3-2R-V
					L	●	128-U3-8L-V	128-U3-2L-V
Cuspid Hk (with Hook)	-7°	8°	0°	3.4 mm	R	●	128-U3-8RK-V	128-U3-2RK-V
					L	●	128-U3-8LK-V	128-U3-2LK-V
Cuspid (without Hook)	0°	8°	0°	3.4 mm	R	●	128-U3-8R0-V	128-U3-2R0-V
					L	●	128-U3-8L0-V	128-U3-2L0-V
Cuspid Hk (with Hook)	0°	8°	0°	3.4 mm	R	●	128-U3-8RK0-V	128-U3-2RK0-V
					L	●	128-U3-8LK0-V	128-U3-2LK0-V
1 st & 2 nd Bicuspids	-7°	0°	0°	3.4 mm	U	●●	128-U4-8U-V	128-U4-2U-V
1 st & 2 nd Bicuspids Hk	-7°	0°	0°	3.4 mm	R	●	128-U4-8RK-V	128-U4-2RK-V
					L	●	128-U4-8LK-V	128-U4-2LK-V

MANDIBULAR MBT² Rx

Tooth	Torque	Angle	Rotation	M/D	R/L	ID	.018	.022
128-L1-2U-V	Anterior	-6°	0°	0°	2.6 mm	U	●●	128-L1-8U-V
Cuspid (without Hook)	-6°	3°	0°	3.4 mm	R	●	128-L3-8R-V	128-L3-2R-V
					L	●	128-L3-8L-V	128-L3-2L-V
Cuspid Hk (with Hook)	-6°	3°	0°	3.4 mm	R	●	128-L3-8RK-V	128-L3-2RK-V
					L	●	128-L3-8LK-V	128-L3-2LK-V
Cuspid (without Hook)	0°	3°	0°	3.4 mm	R	●	128-L3-8R0-V	128-L3-2R0-V
					L	●	128-L3-8L0-V	128-L3-2L0-V
Cuspid Hk (with Hook)	0°	3°	0°	3.4 mm	R	●	128-L3-8RK0-V	128-L3-2RK0-V
					L	●	128-L3-8LK0-V	128-L3-2LK0-V
1 st Bicuspids	-12°	0°	0°	3.4 mm	U	●●	128-L4-8U-V	128-L4-2U-V
1 st Bicuspids Hk	-12°	0°	0°	3.4 mm	R	●	128-L4-8RK-V	128-L4-2RK-V
					L	●	128-L4-8LK-V	128-L4-2LK-V
2 nd Bicuspids	-17°	0°	0°	3.4 mm	U	●●	128-L5-8U-V	128-L5-2U-V
2 nd Bicuspids Hk	-17°	0°	0°	3.4 mm	R	●	128-L5-8RK-V	128-L5-2RK-V
					L	●	128-L5-8LK-V 1	28-L5-2LK-V

BRACKET KITS MBT² Rx

20 brackets per kit	.018	.022
Upper/Lower 5x5 Hook on 3	21K-3-18-01	21K-3-22-011
Upper/Lower 5x5 Hook on 3, 4, & 5	21K-345-18-01	21K-345-22-01

MBT Rx kits include 0° torque on the cuspids - custom kits available upon request.

Brackets are packaged with 1 case oriented in a kit box.

Brackets are also available in kits of 5 or 10 cases oriented in a carousel dispenser.



Agility™ Twin Bracket
kit, 20 brackets
oriented in a kit box

