




## Springs & Auxiliaries/ Ressorts & Auxiliaires

# INITIALLOY Springs



## INITIALLOY Stop Wound™ Open Springs

| Description | Force | Color  | Internal bore | P/N          |
|-------------|-------|--------|---------------|--------------|
| Extra-Light | 50 g  | Black  | .036"         | 60-0112-0000 |
| Light       | 100 g | Blue   | .036"         | 60-0113-0000 |
| Medium      | 150 g | Yellow | .036"         | 60-0114-0000 |
| Medium      | 150 g | Yellow | .045"         | 60-0114-0045 |
| Heavy       | 200 g | Red    | .036"         | 60-0115-0000 |

3 stuks/verpakking - 3 pcs/ emballage



## INITIALLOY Open Springs

| Description | Force | Color  | P/N          |
|-------------|-------|--------|--------------|
| Medium      | 150 g | Yellow | 60-0124-0000 |

10 stuks/verpakking - 10 pcs/emballage



## NITI Open Spring (in Spool)

|                                                                    | P/N          |
|--------------------------------------------------------------------|--------------|
| INITIALLOY Open Spring in Spool, Medium Force (150 g), .010"x.036" | 60-0144-0000 |
| Standard Open Spring in Spool .010" x .030"                        | 60-0140-2013 |



## INITIALLOY Closed Springs

| Description         | Force | Color  | P/N          |
|---------------------|-------|--------|--------------|
| Light               | 100 g | Blue   | 60-0223-0000 |
| Medium              | 150 g | Yellow | 60-0224-0000 |
| Heavy               | 200 g | Red    | 60-0225-0000 |
| Incisors Retraction | 300 g | White  | 60-0227-0000 |

10 stuks/verpakking - 10 pcs/emballage

# STAINLESS STEEL Springs



## OPEN SPRING (91cm / 3' Spool)

| Dimension   | P/N          |
|-------------|--------------|
| .008"x.022" | 61-0140-0082 |
| .010"x.030" | 61-0140-0013 |
| .010"x.040" | 61-0140-0014 |

## CLOSED SPRING (91cm / 3' Spool)

| Dimension   | P/N          |
|-------------|--------------|
| .008"x.022" | 61-0240-0082 |
| .010"x.030" | 61-0240-0013 |



# Constant Force NiTi Springs

| Item                                                                                                      | Force | ID      | Wire Size | 9mm     | 12mm    |
|-----------------------------------------------------------------------------------------------------------|-------|---------|-----------|---------|---------|
|  <b>Extra Light (XL)</b> | 50g   | 1 Hole  | .007x.038 | 226-290 | 226-295 |
|  <b>Light (L)</b>        | 100g  | 2 Holes | .009x.038 | 226-291 | 226-296 |
|  <b>Medium (M)</b>       | 150g  | 3 Holes | .010x.038 | 226-292 | 226-297 |
|  <b>Heavy (H)</b>        | 200g  | 4 Holes | .011x.042 | 226-293 | 226-298 |
|  <b>Extra Heavy (XH)</b> | 250g  | 5 Holes | .012x.042 | 226-294 | 226-299 |

Constant Force Spring Kit:

\*Full kit includes 1 each of: 32mm, 34mm, 36mm, 38mm, and 40mm 5 of each size, 9 mm & 1 2 mm (50 total) 226-289

5 stuks/verpakking - 5 pcs/emballage



| (Passive Spring Length) |      | (Force Rating) |               | Extension Measurement Landmarks (mm) / Subsequent Force Values (gms) |      |      |      |      |      |      |      |      |
|-------------------------|------|----------------|---------------|----------------------------------------------------------------------|------|------|------|------|------|------|------|------|
| Length                  | Size | Force          | Order Numbers | 13mm                                                                 | 15mm | 17mm | 19mm | 22mm | 27mm | 31mm | 36mm | 41mm |
| Extra Short<br>(11mm)   | A1   | Light (150g)   | CCOT45LT      | 95                                                                   | 125  | 150  | -    | -    | -    | -    | -    | -    |
|                         | A2   | Medium (200g)  | CCOT45MD      | 130                                                                  | 170  | 200  | -    | -    | -    | -    | -    | -    |
|                         | A3   | Heavy (250g)   | CCOT45HV      | 150                                                                  | 210  | 250  | -    | -    | -    | -    | -    | -    |
| Short<br>(13mm)         | B1   | Light (150g)   | CCOT6LT       | -                                                                    | 95   | 120  | 130  | 150  | -    | -    | -    | -    |
|                         | B3   | Medium (200g)  | CCOT6MD       | -                                                                    | 125  | 150  | 170  | 200  | -    | -    | -    | -    |
|                         | B3   | Heavy (250g)   | CCOT6HV       | -                                                                    | 150  | 180  | 200  | 250  | -    | -    | -    | -    |
| Medium<br>(15mm)        | C1   | Light (150g)   | CCOT9LT       | -                                                                    | -    | 85   | 100  | 110  | 130  | 150  | -    | -    |
|                         | C2   | Medium (200g)  | CCOT9MD       | -                                                                    | -    | 115  | 140  | 150  | 180  | 200  | -    | -    |
|                         | C3   | Heavy (250g)   | CCOT9HV       | -                                                                    | -    | 145  | 165  | 175  | 220  | 250  | -    | -    |
| Long<br>(17mm)          | D1   | Light (150g)   | CCOT12LT      | -                                                                    | -    | -    | 85   | 95   | 110  | 120  | 130  | 150  |
|                         | D2   | Medium (200g)  | CCOT12MD      | -                                                                    | -    | -    | 110  | 125  | 150  | 170  | 190  | 200  |
|                         | D3   | Heavy (250g)   | CCOT12HV      | -                                                                    | -    | -    | 130  | 150  | 175  | 200  | 230  | 250  |

10 stuks/verpakking - 10 pcs/emballage



## G4™ NiTi Closed Coil Springs

|             | Force 21" | Spool     |
|-------------|-----------|-----------|
| Ultra Light | 100 gm    | CCNT21075 |
| Light       | 150 gm    | CCNT21150 |
| Medium      | 200 gm    | CCNT21200 |
| Heavy       | 250 gm    | CCNT21250 |



## G4™ NiTi Open Coil Springs

| Wire size | I.D. | Force     | Force when compressed 25% (.250")<br>(Force in grams) | Force when compressed 50% (.500")<br>(Force in grams) | 7" Straight 21" | Spooled    |
|-----------|------|-----------|-------------------------------------------------------|-------------------------------------------------------|-----------------|------------|
| .010"     | .030 | Light     | 120                                                   | 260                                                   | OCNT10307       | OCNT103021 |
| .010"     | .036 | Med-Light | 90                                                    | 200                                                   | OCNT10367       | OCNT103621 |
| .012"     | .030 | Medium    | 220                                                   | 480                                                   | OCNT12307       | OCNT123021 |

3 stuks/verpakking - 3 pcs/emballage



## G4™ NiTi Open Coil Springs

| I.D. | Force |        |           |
|------|-------|--------|-----------|
| .036 | 100gm | Light  | OCNT15100 |
| .036 | 150gm | Medium | OCNT15150 |
| .036 | 200gm | Heavy  | OCNT15200 |

10 stuks/verpakking - 10 pcs/emballage

**Springs/Auxiliaires**

# G4 NiTi Open Coil Springs



## G4 Niti Molar Distalizing Open Coil Springs

|            |                 |
|------------|-----------------|
| 7" Lengths | OCNT7045 (3pcs) |
| 21" Spool  | OCNT2145        |

# Stainless Steel Coil Springs



## STAINLESS STEEL Open Coil Springs (spool)

|           |          |
|-----------|----------|
| .010x.028 | SSOC1028 |
| .010x.030 | SSOC1030 |
| .010x.036 | SSOC1036 |



## STAINLESS STEEL Closed Coil (spool)

|           |          |
|-----------|----------|
| .010x.030 | SSCC1030 |
|-----------|----------|

# BONDABLE Lingual Attachments

## BONDABLE Buttons

P/N



Flat base

42-0020-0000



Curved base

42-0010-0000

## BONDABLE Cleats

P/N



Flat base

42-0060-0000



Curved base

42-0050-0000



Eyelet

42-0070-0000

## SURGICAL Buttons

P/N



Curved base

42-0030-0000



Flat base

42-0040-0000



10 stuks/verpakking - 10 pcs/emballage

## V-SLOT OMNI HOOK

| Description                 | P/N          |
|-----------------------------|--------------|
| Hook for V-Slot (Bag of 50) | 69-1410-0000 |

10 stuks/verpakking - 10 pcs/emballage

# BONDABLE Lingual Attachments

## COIL EYELETS

**M006B-01** - 10 pcs



## CERAMIC BUTTONS

**FCB-00-01** - 10 pcs



## CERAMIC HOOKS

**FCB-00-01HK** - 10 pcs



## COMPOSITE BUTTONS

**270-204** - 10 pcs



## MIM DIRECT BOND EYELET

**270-203** - 10 pcs



## GOLD ERUPTION CHAIN

**42-3150-0000**

Gold Eruption Chain  
round base - 5 pcs



**42-3160-0000**

Gold Eruption Chain  
rectangular base - 5 pcs



## TONGUE TRAINER

**270-205**


Tongue Trainer - 10 pcs




## BITE BUMPER

**BB-00-01**


BITE BUMPERS - aesthetic- 10 pcs



1. Put the placement-jig against the incisal edge



2. Remove the jig



3. Bite-Bumper ready

## MINI-MOLD SYSTEM

**700-100M+** Starterkit

Kit includes: 24 Mini-Mold tips (3 of each style),  
2 tip handles, detailed step-by-step instructions.



**Bite Ramp**

**700-117** 3mm  
8 tips, 1 handle  
**700-117L** 5mm  
8 tips, 1 handle



**Button**

**700-004**  
10 tips,  
1 handle



**Habit Control**

**700-116**  
10 tips,  
1 handle




**Molar Pad**

**700-007**  
10 tips,  
1 handle

## GOODMAN TORQUING SPRING™

- Torque veertje
- .013 roestvrij staal draad
- 125 tot 150 gram labiale of linguale root torque
- Crimpt op elke rechthoekige draad van .016x.022 tot .019x.025
- Verkrijgbaar in twee maten: Small & Large
- De small-versie is voor de laterale tanden upper/lower en de onderincisiven. De large voor al de andere tanden (niet voor de molaren)
- Verpakt per 5 stuks



- Ressort de torsion
- fil d'acier inoxydable .013
- 125 à 150 grammes de torque labial ou lingual
- Sertissage pour tout fil rectangulaire de .016x.022 à .019x.025
- Disponible en deux tailles : Small & Large
- La version small est destinée aux dents latérales supérieures/inférieures et aux incisives inférieures. La version large pour toutes les autres dents (pas pour les molaires).
- Emballé par 5 pièces



### GTS Complete kit 270-104

- 5 pcs small
- 5 pcs large
- 1 crimping plier

- 270-100**  
5 pcs small
- 270-101**  
5 pcs large



GTS Crimping plier  
800-087

## D2 APPLIANCE 271-600 one pair

- One size fits all om uw voorraad sterk te verminderen
- Laag profielontwerp voor meer comfort voor de patiënt
- Voor bilaterale of unilaterale behandeling
- Voor gemengd en permanent gebit
- Eenvoudig te gebruiken en gemakkelijk te plaatsen
- De totale behandelingstijd verkorten
- Gepatenteerd voor zowel klasse II als klasse III behandeling

- Une taille unique pour réduire considérablement votre inventaire
- Un profil bas pour le confort du patient
- Pour un traitement bilatéral ou unilatéral
- Pour les dents mixtes et permanentes
- Simple à utiliser et facile à placer
- Réduire la durée globale du traitement
- Breveté pour un traitement de classe II et de classe III



## D2 APPLIANCE CLEAR 271-600C





- Simple
- Precise
- Adaptable
- Efficient

Youtube introduction by  
Dr.Yoav Mazor




# MULTI-FA

Multi functional Appliance  
for CL II and CL III Correction

## Multi FA clinical uses:

1. CL II, CL III malocclusion correction
2. Open / regain spaces
3. Closing spaces



Double hinge for an  
anatomical positionning  
without compromise



Improvement of  
the biomechanical control



Fix the hinge upon  
treatment plan  
requirements



**One appliance  
for 90% cases!**

## MULTI-FA for tubes, adaptable on patients

Can be used with open coil



mesio – distal stop made with a flow



or a crimpable stop

**92-5001-0000**  
MULTI-FA Regular 18-28mm  
(left & right) - 2 pcs

**92-5002-0000**  
MULTI FA Short 12-17mm  
(Left & Right) - 2 pcs

**92-5003-0000**  
MULTI FA Long 23-36mm  
(Left & Right) - 2 pcs

**92-5004-0000**  
MULTI FA Band-R  
(Left & Right) - 2 pcs