




## Springs & Auxiliaries/ Ressorts & Auxiliaires

# INITIALLOY Springs



## INITIALLOY Stop Wound™ Open Springs

Description	Force	Color	Internal bore	P/N
Extra-Light	50 g	Black	.036"	60-0112-0000
Light	100 g	Blue	.036"	60-0113-0000
Medium	150 g	Yellow	.036"	60-0114-0000
Medium	150 g	Yellow	.045"	60-0114-0045
Heavy	200 g	Red	.036"	60-0115-0000

3 stuks/verpakking - 3 pcs/ emballage



## INITIALLOY Open Springs

Description	Force	Color	P/N
Medium	150 g	Yellow	60-0124-0000

10 stuks/verpakking - 10 pcs/emballage



## NITI Open Spring (in Spool)

	P/N
INITIALLOY Open Spring in Spool, Medium Force (150 g), .010"x.036"	60-0144-0000
Standard Open Spring in Spool .010" x .030"	60-0140-2013



## INITIALLOY Closed Springs

Description	Force	Color	P/N
Light	100 g	Blue	60-0223-0000
Medium	150 g	Yellow	60-0224-0000
Heavy	200 g	Red	60-0225-0000
Incisors Retraction	300 g	White	60-0227-0000

10 stuks/verpakking - 10 pcs/emballage

# STAINLESS STEEL Springs



## OPEN SPRING (91cm / 3' Spool)






Dimension	P/N
.008"x.022"	61-0140-0082
.010"x.030"	61-0140-0013
.010"x.040"	61-0140-0014

## CLOSED SPRING (91cm / 3' Spool)

Dimension	P/N
.008"x.022"	61-0240-0082
.010"x.030"	61-0240-0013



# Constant Force NiTi Springs

Item	Force	ID	Wire Size	9mm	12mm
 <b>Extra Light (XL)</b>	50g	1 Hole	.007x.038	226-290	226-295
 <b>Light (L)</b>	100g	2 Holes	.009x.038	226-291	226-296
 <b>Medium (M)</b>	150g	3 Holes	.010x.038	226-292	226-297
 <b>Heavy (H)</b>	200g	4 Holes	.011x.042	226-293	226-298
 <b>Extra Heavy (XH)</b>	250g	5 Holes	.012x.042	226-294	226-299

Constant Force Spring Kit:

\*Full kit includes 1 each of: 32mm, 34mm, 36mm, 38mm, and 40mm 5 of each size, 9 mm & 12 mm (50 total) 226-289

5 stuks/verpakking - 5 pcs/emballage



(Passive Spring Length)	(Force Rating)	Extension Measurement Landmarks (mm) / Subsequent Force Values (gms)											
		Length	Size	Force	Order Numbers	13mm	15mm	17mm	19mm	22mm	27mm	31mm	36mm
<b>Extra Short</b> (11mm)	<b>A1</b>	Light (150g)	CCOT45LT	95	125	150	-	-	-	-	-	-	-
	<b>A2</b>	Medium (200g)	CCOT45MD	130	170	200	-	-	-	-	-	-	
	<b>A3</b>	Heavy (250g)	CCOT45HV	150	210	250	-	-	-	-	-	-	
<b>Short</b> (13mm)	<b>B1</b>	Light (150g)	CCOT6LT	-	95	120	130	150	-	-	-	-	
	<b>B3</b>	Medium (200g)	CCOT6MD	-	125	150	170	200	-	-	-	-	
	<b>B3</b>	Heavy (250g)	CCOT6HV	-	150	180	200	250	-	-	-	-	
<b>Medium</b> (15mm)	<b>C1</b>	Light (150g)	CCOT9LT	-	-	85	100	110	130	150	-	-	
	<b>C2</b>	Medium (200g)	CCOT9MD	-	-	115	140	150	180	200	-	-	
	<b>C3</b>	Heavy (250g)	CCOT9HV	-	-	145	165	175	220	250	-	-	
<b>Long</b> (17mm)	<b>D1</b>	Light (150g)	CCOT12LT	-	-	-	85	95	110	120	130	150	
	<b>D2</b>	Medium (200g)	CCOT12MD	-	-	-	110	125	150	170	190	200	
	<b>D3</b>	Heavy (250g)	CCOT12HV	-	-	-	130	150	175	200	230	250	

10 stuks/verpakking - 10 pcs/emballage



## G4™ NiTi Closed Coil Springs

	Force 21"	Spool
Ultra Light	100 gm	CCNT21075
Light	150 gm	CCNT21150
Medium	200 gm	CCNT21200
Heavy	250 gm	CCNT21250



## G4™ NiTi Open Coil Springs

Wire size	I.D.	Force	Force when compressed 25% (.250") (Force in grams)	Force when compressed 50% (.500") (Force in grams)	7" Straight 21"	Spooled
.010"	.030	Light	120	260	OCNT10307	OCNT103021
.010"	.036	Med-Light	90	200	OCNT10367	OCNT103621
.012"	.030	Medium	220	480	OCNT12307	OCNT123021

3 stuks/verpakking - 3 pcs/emballage




## G4™ NiTi Open Coil Springs

I.D.	Force		
.036	100gm	Light	OCNT15100
.036	150gm	Medium	OCNT15150
.036	200gm	Heavy	OCNT15200



10 stuks/verpakking - 10 pcs/emballage

Springs/Auxiliaires

# G4 NiTi Open Coil Springs

	<b>G4 Niti Molar Distalizing Open Coil Springs</b>	
	7" Lengths	OCNT7045 (3pcs)
	21" Spool	OCNT2145

# Stainless Steel Coil Springs

	<b>STAINLESS STEEL Open Coil Springs (spool)</b>	
	.010x.028	SSOC1028
	.010x.030	SSOC1030
	.010x.036	SSOC1036
	<b>STAINLESS STEEL Closed Coil (spool)</b>	
	.010x.030	S SCC1030

# BONDABLE Lingual Attachments

<b>BONDABLE Buttons</b>		P/N
		Flat base 42-0020-0000
		Curved base 42-0010-0000
<b>BONDABLE Cleats</b>		P/N
		Flat base 42-0060-0000
		Curved base 42-0050-0000
		Eyelet 42-0070-0000
<b>SURGICAL Buttons</b>		P/N
	 1.3	Curved base 42-0030-0000
	 1.3	Flat base 42-0040-0000
10 stuks/verpakking - 10 pcs/emballage		
	<b>V-SLOT OMNI HOOK</b>	
	Description	P/N
	Hook for V-Slot (Bag of 50)	69-1410-0000

10 stuks/verpakking - 10 pcs/emballage

# BONDABLE Lingual Attachments

## COIL EYELETS

**M006B-01** - 10 pcs



## CERAMIC BUTTONS

**FCB-00-01** - 10 pcs



## CERAMIC HOOKS

**FCB-00-01HK** - 10 pcs



## COMPOSITE BUTTONS

**270-204** - 10 pcs



## MIM DIRECT BOND EYELET

**270-203** - 10 pcs



## GOLD ERUPTION CHAIN

**42-3150-0000**

Gold Eruption Chain  
round base - 5 pcs



**42-3160-0000**

Gold Eruption Chain  
rectangular base - 5 pcs



## TONGUE TRAINER

**270-205**


Tongue Trainer - 10 pcs




## BITE BUMPER

**BB-00-01**


BITE BUMPERS - aesthetic- 10 pcs



1. Put the placement-jig against the incisal edge



2. Remove the jig



3. Bite-Bumper ready

## MINI-MOLD SYSTEM

**700-100M+** Starterkit

Kit includes: 24 Mini-Mold tips (3 of each style),  
2 tip handles, detailed step-by-step instructions.



**Bite Ramp**

**700-117** 3mm  
8 tips, 1 handle  
**700-117L** 5mm  
8 tips, 1 handle



**Button**

**700-004**  
10 tips,  
1 handle



**Habit Control**

**700-116**  
10 tips,  
1 handle



**Molar Pad**

**700-007**  
10 tips,  
1 handle

## GOODMAN TORQUING SPRING™

- Torque veertje
- .013 roestvrij staal draad
- 125 tot 150 gram labiale of linguale root torque
- Crimpt op elke rechthoekige draad van .016x.022 tot .019x.025
- Verkrijgbaar in twee maten: Small & Large
- De small-versie is voor de laterale tanden upper/lower en de onderincisiven. De large voor al de andere tanden. (niet voor de molaren)
- Verpakt per 5 stuks



- Ressort de torsion
- fil d'acier inoxydable .013
- 125 à 150 grammes de torque labial ou lingual
- Sertissage pour tout fil rectangulaire de .016x.022 à .019x.025
- Disponible en deux tailles : Small & Large
- La version small est destinée aux dents latérales supérieures/inférieures et aux incisives inférieures. La version large pour toutes les autres dents (pas pour les molaires).
- Emballé par 5 pièces



### GTS Complete kit 270-104

- 5 pcs small
- 5 pcs large
- 1 crimping plier

- 270-100**  
5 pcs small
- 270-101**  
5 pcs large



GTS Crimping plier  
800-087

## D2 APPLIANCE 271-600 one pair

- One size fits all om uw voorraad sterk te verminderen
- Laag profielontwerp voor meer comfort voor de patiënt
- Voor bilaterale of unilaterale behandeling
- Voor gemengd en permanent gebit
- Eenvoudig te gebruiken en gemakkelijk te plaatsen
- De totale behandelingstijd verkorten
- Gepatenteerd voor zowel klasse II als klasse III behandeling

- Une taille unique pour réduire considérablement votre inventaire
- Un profil bas pour le confort du patient
- Pour un traitement bilatéral ou unilatéral
- Pour les dents mixtes et permanentes
- Simple à utiliser et facile à placer
- Réduire la durée globale du traitement
- Breveté pour un traitement de classe II et de classe III

ONE SIZE  
FITS ALL



## D2 APPLIANCE CLEAR 271-600C







- Simple
- Precise
- Adaptable
- Efficient

Youtube introduction by  
Dr. Yoav Mazor




# MULTI-FA

Multi functional Appliance  
for CL II and CL III Correction

### Multi FA clinical uses:

1. CL II, CL III malocclusion correction
2. Open / regain spaces
3. Closing spaces




Double hinge for an  
anatomical positioning  
without compromise



Improvement of  
the biomechanical control




Fix the hinge upon  
treatment plan  
requirements



One appliance  
for 90% cases!

### MULTI-FA for tubes, adaptable on patients

Can be used with open coil



mesio – distal stop made with a flow



or a crimpable stop

**92-5001-0000**  
MULTI-FA Regular 18-28mm  
(left & right) - 2 pcs

**92-5002-0000**  
MULTI FA Short 12-17mm  
(Left & Right) - 2 pcs

**92-5003-0000**  
MULTI FA Long 23-36mm  
(Left & Right) - 2 pcs

**92-5004-0000**  
MULTI FA Band-R  
(Left & Right) - 2 pcs