




Springs & Auxiliaries/ Ressorts & Auxiliaires

INITIALLOY Springs



INITIALLOY Stop Wound™ Open Springs

Description	Force	Color	Internal bore	P/N
Extra-Light	50 g	Black	.036"	60-0112-0000
Light	100 g	Blue	.036"	60-0113-0000
Medium	150 g	Yellow	.036"	60-0114-0000
Medium	150 g	Yellow	.045"	60-0114-0045
Heavy	200 g	Red	.036"	60-0115-0000

3 stuks/verpakking - 3 pcs/ emballage



INITIALLOY Open Springs

Description	Force	Color	P/N
Medium	150 g	Yellow	60-0124-0000

10 stuks/verpakking - 10 pcs/emballage



NITI Open Spring (in Spool)

	P/N
INITIALLOY Open Spring in Spool, Medium Force (150 g), .010"x.036"	60-0144-0000
Standard Open Spring in Spool .010" x .030"	60-0140-0013



INITIALLOY Closed Springs

Description	Force	Color	P/N
Light	100 g	Blue	60-0223-0000
Medium	150 g	Yellow	60-0224-0000
Heavy	200 g	Red	60-0225-0000
Incisors Retraction	300 g	White	60-0227-0000

10 stuks/verpakking - 10 pcs/emballage

STAINLESS STEEL Springs



OPEN SPRING (91cm / 3' Spool)






Dimension	P/N
.008"x.022"	61-0140-0082
.010"x.030"	61-0140-0013
.010"x.040"	61-0140-0014

CLOSED SPRING (91cm / 3' Spool)

Dimension	P/N
.008"x.022"	61-0240-0082
.010"x.030"	61-0240-0013



Constant Force NiTi Springs

Item	Force	ID	Wire Size	9mm	12mm
 Extra Light (XL)	50g	1 Hole	.007x.038	226-290	226-295
 Light (L)	100g	2 Holes	.009x.038	226-291	226-296
 Medium (M)	150g	3 Holes	.010x.038	226-292	226-297
 Heavy (H)	200g	4 Holes	.011x.042	226-293	226-298
 Extra Heavy (XH)	250g	5 Holes	.012x.042	226-294	226-299

Constant Force Spring Kit:

*Full kit includes 1 each of: 32mm, 34mm, 36mm, 38mm, and 40mm 5 of each size, 9 mm & 12 mm (50 total) 226-289

5 stuks/verpakking - 5 pcs/emballage



(Passive Spring Length)	(Force Rating)	Extension Measurement Landmarks (mm) / Subsequent Force Values (gms)											
		Length	Size	Force	Order Numbers	13mm	15mm	17mm	19mm	22mm	27mm	31mm	36mm
Extra Short (11mm)	A1	Light (150g)	CCOT45LT	95	125	150	-	-	-	-	-	-	-
	A2	Medium (200g)	CCOT45MD	130	170	200	-	-	-	-	-	-	-
	A3	Heavy (250g)	CCOT45HV	150	210	250	-	-	-	-	-	-	-
Short (13mm)	B1	Light (150g)	CCOT6LT	-	95	120	130	150	-	-	-	-	-
	B3	Medium (200g)	CCOT6MD	-	125	150	170	200	-	-	-	-	-
	B3	Heavy (250g)	CCOT6HV	-	150	180	200	250	-	-	-	-	-
Medium (15mm)	C1	Light (150g)	CCOT9LT	-	-	85	100	110	130	150	-	-	-
	C2	Medium (200g)	CCOT9MD	-	-	115	140	150	180	200	-	-	-
	C3	Heavy (250g)	CCOT9HV	-	-	145	165	175	220	250	-	-	-
Long (17mm)	D1	Light (150g)	CCOT12LT	-	-	-	85	95	110	120	130	150	-
	D2	Medium (200g)	CCOT12MD	-	-	-	110	125	150	170	190	200	-
	D3	Heavy (250g)	CCOT12HV	-	-	-	130	150	175	200	230	250	-

10 stuks/verpakking - 10 pcs/emballage



G4™ NiTi Closed Coil Springs

	Force 21"	Spool
Ultra Light	100 gm	CCNT21075
Light	150 gm	CCNT21150
Medium	200 gm	CCNT21200
Heavy	250 gm	CCNT21250



G4™ NiTi Open Coil Springs

Wire size	I.D.	Force	Force when compressed 25% (.250") (Force in grams)	Force when compressed 50% (.500") (Force in grams)	7" Straight 21"	Spooled
.010"	.030	Light	120	260	OCNT10307	OCNT103021
.010"	.036	Med-Light	90	200	OCNT10367	OCNT103621
.012"	.030	Medium	220	480	OCNT12307	OCNT123021

3 stuks/verpakking - 3 pcs/emballage




G4™ NiTi Open Coil Springs

I.D.	Force		
.036	100gm	Light	OCNT15100
.036	150gm	Medium	OCNT15150
.036	200gm	Heavy	OCNT15200



10 stuks/verpakking - 10 pcs/emballage

Springs/Auxiliaires

G4 NiTi Open Coil Springs

	G4 Niti Molar Distalizing Open Coil Springs	
	7" Lengths	OCNT7045 (3pcs)
	21" Spool	OCNT2145

Stainless Steel Coil Springs

	STAINLESS STEEL Open Coil Springs (spool)	
	.010x.028	SSOC1028
	.010x.030	SSOC1030
	.010x.036	SSOC1036
	STAINLESS STEEL Closed Coil (spool)	
	.010x.030	S SCC1030

BONDABLE Lingual Attachments

BONDABLE Buttons		P/N
		Flat base 42-0020-0000
		Curved base 42-0010-0000
BONDABLE Cleats		P/N
		Flat base 42-0060-0000
		Curved base 42-0050-0000
		Eyelet 42-0070-0000
SURGICAL Buttons		P/N
	 1.3	Curved base 42-0030-0000
	 1.3	Flat base 42-0040-0000
10 stuks/verpakking - 10 pcs/emballage		
	V-SLOT OMNI HOOK	
	Description	P/N
	Hook for V-Slot (Bag of 50)	69-1410-0000

10 stuks/verpakking - 10 pcs/emballage

BONDABLE Lingual Attachments

COIL EYELETS

M006B-01 - 10 pcs



CERAMIC BUTTONS

FCB-00-01 - 10 pcs



CERAMIC HOOKS

FCB-00-01HK - 10 pcs



COMPOSITE BUTTONS

270-204 - 10 pcs



MIM DIRECT BOND EYELET

270-203 - 10 pcs



GOLD ERUPTION CHAIN

42-3150-0000

Gold Eruption Chain
round base - 5 pcs



42-3160-0000

Gold Eruption Chain
rectangular base - 5 pcs



TONGUE TRAINER

270-205


Tongue Trainer - 10 pcs




BITE BUMPER

BB-00-01


BITE BUMPERS - aesthetic- 10 pcs



1. Put the placement-jig against the incisal edge



2. Remove the jig



3. Bite-Bumper ready

MINI-MOLD SYSTEM

700-100M+ Starterkit

Kit includes: 24 Mini-Mold tips (3 of each style),
2 tip handles, detailed step-by-step instructions.



Bite Ramp

700-117 3mm
8 tips, 1 handle
700-117L 5mm
8 tips, 1 handle



Button

700-004
10 tips,
1 handle



Habit Control

700-116
10 tips,
1 handle



Molar Pad

700-007
10 tips,
1 handle

GOODMAN TORQUING SPRING™

- Torque veertje
- .013 roestvrij staal draad
- 125 tot 150 gram labiale of linguale root torque
- Crimpt op elke rechthoekige draad van .016x.022 tot .019x.025
- Verkrijgbaar in twee maten: Small & Large
- De small-versie is voor de laterale tanden upper/lower en de onderincisiven. De large voor al de andere tanden. (niet voor de molaren)
- Verpakt per 5 stuks



- Ressort de torsion
- fil d'acier inoxydable .013
- 125 à 150 grammes de torque labial ou lingual
- Sertissage pour tout fil rectangulaire de .016x.022 à .019x.025
- Disponible en deux tailles : Small & Large
- La version small est destinée aux dents latérales supérieures/inférieures et aux incisives inférieures. La version large pour toutes les autres dents (pas pour les molaires).
- Emballé par 5 pièces



GTS Complete kit 270-104

- 5 pcs small
- 5 pcs large
- 1 crimping plier

- 270-100**
5 pcs small
- 270-101**
5 pcs large



GTS Crimping plier
800-087

D2 APPLIANCE 271-600 one pair

- One size fits all om uw voorraad sterk te verminderen
- Laag profielontwerp voor meer comfort voor de patiënt
- Voor bilaterale of unilaterale behandeling
- Voor gemengd en permanent gebit
- Eenvoudig te gebruiken en gemakkelijk te plaatsen
- De totale behandelingstijd verkorten
- Gepatenteerd voor zowel klasse II als klasse III behandeling

- Une taille unique pour réduire considérablement votre inventaire
- Un profil bas pour le confort du patient
- Pour un traitement bilatéral ou unilatéral
- Pour les dents mixtes et permanentes
- Simple à utiliser et facile à placer
- Réduire la durée globale du traitement
- Breveté pour un traitement de classe II et de classe III

ONE SIZE
FITS ALL



D2 APPLIANCE CLEAR 271-600C





- Simple
- Precise
- Adaptable
- Efficient

Youtube introduction by
Dr. Yoav Mazor




MULTI-FA

Multi functional Appliance
for CL II and CL III Correction

Multi FA clinical uses:

1. CL II, CL III malocclusion correction
2. Open / regain spaces
3. Closing spaces



Double hinge for an
anatomical positioning
without compromise



Improvement of
the biomechanical control



Fix the hinge upon
treatment plan
requirements



**One appliance
for 90% cases!**

MULTI-FA for tubes, adaptable on patients

Can be used with open coil



mesio – distal stop made with a flow



or a crimpable stop

92-5001-0000
MULTI-FA Regular 18-28mm
(left & right) - 2 pcs

92-5002-0000
MULTI FA Short 12-17mm
(Left & Right) - 2 pcs

92-5003-0000
MULTI FA Long 23-36mm
(Left & Right) - 2 pcs

92-5004-0000
MULTI FA Band-R
(Left & Right) - 2 pcs